

Fall Prevention Systems



PROFESSIONAL MEDICAL, INC.

Fall Prevention Alarms, Sensor Pads and Systems

STRING ALARM



The Compliance String Alarm is our most economical alarm offering. It is rich in features and is made of durable materials that can withstand the rigors of institutional use.

- Monitor is always on, making it impossible for staff to forget to turn it on
- Unique alarm tone provides differentiation from other alarms in your facility
- Low battery indicator is located on the front of the monitor for easy visibility
- Nylon strap easily hangs from the push handle of a wheelchair
- Optional stainless steel mounting clip for mounting to wheelchair back upholstery or storage pockets
- Available with optional protective boot to protect the monitor from accidental drops (600-ALPB)

Item #
600-AA1

MOBILE ALARM



The feature-rich, low-cost Compliance Mobile Alarm provides unparalleled value and is an excellent entry-level upgrade from string monitors.

- Monitor is always on, making it impossible for staff to forget to turn it on
- Alarm automatically resets when pressure is re-applied to the sensor pad
- Low battery indicator is located on the front of the monitor for easy visibility
- Strain relief tab helps reduce strain on sensor pad cable and protects the connector
- Snap-to-close battery cover also has a screw for increased tamper resistance
- Comes with an unbreakable stainless steel mounting clip
- Available with optional protective boot to protect the monitor from accidental drops (600-ALPB)
- Compatible with our universal chair and bed sensors, seat belt systems

Item #
600-AL1

FLEX ALARM



The Compliance Flex Alarm is compatible with nurse call systems and other monitoring tools.

- Monitor is always on, making it impossible for staff to forget to turn it on
- Alarm automatically resets when pressure is re-applied to the sensor pad
- Alarm indicator on the front of the monitor flashes when in alarm mode
- All alarm settings (Hi/Lo/remote volume, floor/sensor pad selection) are located inside the battery compartment for tamper resistance
- Strain relief tab helps reduce strain on sensor pad cable and protects the connector
- Unbreakable stainless steel mounting clip
- Optional AC adapter available
- Available with optional protective boot to protect the monitor from accidental drops (600-ALPB)
- Compatible with our universal chair and bed sensors, seat belt systems and floor sensor pad (600-ALFP-365)

Item #
600-AL2

ULTRA ALARM



The Compliance 3-in-1 Ultra Alarm combines many of the great features of our other alarms with monitoring tool compatibility for ultimate flexibility and convenience.

- Monitor is always on, making it impossible for staff to forget to turn it on
- Alarm automatically resets when pressure is re-applied to the sensor pad
- All alarm settings are located inside the battery compartment (Hi/Lo/remote volume, 0-2 second delay)
- Snap-to-close battery cover also has a screw for increased tamper resistance
- Unbreakable stainless steel mounting clip
- Nurse call compatible (also includes remote alarm option)
- Includes protective boot to protect the monitor from accidental drops
- Compatible with our universal chair and bed sensors, seat belt systems and floor sensor pad (600-ALFP-365)

Item #
600-AL3

UNIVERSAL ONE-YEAR SENSOR PADS



Moisture Resistant

One-year moisture-resistant sensor pads for chairs and beds work with many brands of fall monitors.

- Available with one-year moisture-resistant warranty period
- Measures 10" W x 15" L (chair) or 10" W x 28" L (bed)
- 68" sensor cable provides more alarm mounting options (bed only)

Item #
600-ALCP-365
Chair Sensor Pad; 10" W x 15" L
600-ALPB-365
Bed Sensor Pad; 10" W x 28" L

UNIVERSAL 45-DAY SENSOR PADS

Universal 45-day sensor pads for chairs and beds work with many brands of fall monitors.

- Available with 45-day warranty period
- Measures 10" W x 15" L (chair) or 10" W x 28" L (bed)
- 68" sensor cable provides more alarm mounting options (bed only)

Item #
600-ALCP-45
Chair Sensor Pad; 10" W x 15" L
600-ALBP-45
Bed Sensor Pad; 10" W x 28" L

FALL PREVENTION SYSTEMS



Moisture Resistant

Fall prevention systems feature a Mobile Alarm and one-year moisture-resistant chair or bed sensor pad.

Item #
600-AL365-CSYS
Chair System
600-AL365-BSYS
Bed System

FLOOR SENSOR PAD



Weight-sensing floor mats can be used at the bedside or at a doorway to alert caregivers of resident activities.

- Heavy duty and ideal for institutional use
- For use with the Ultra (600-AL3) and Flex (600-AL2) alarms
- 22" x 39"

Item #
600-ALFP-365

SEATBELT MONITORING SYSTEMS



Hook-and-loop and buckle-style seatbelt monitoring systems alarm as soon as a resident attempts to exit his chair. The large red loop on the hook-and-loop seatbelt is easy to release, even for individuals with limited finger dexterity.

- Mounts to most wheelchairs using existing hardware
- Compatible with the Ultra (600-AL3), Flex (600-AL2) and Mobile (600-AL1) alarms

Item #
600-ALBV
Hook-and-Loop Seatbelt Monitoring System
600-ALBB
Buckle Seatbelt Monitoring System



Sensor cords are reinforced to meet the demands of institutional use.

Accessories

BED CLIP



Durable bed clips are designed to hold a fall alarm securely to a bed headboard or footboard.

- Durable, long-lasting design
- Easy to install

Item #

600-ALR-BC

AC ADAPTER



Although the Flex alarm (600-AL2) can operate on batteries alone, this AC adapter saves on battery life when the alarm is used with a bed system.

- Durable, long-lasting design
- 8 feet in length
- Easily reaches power sources

Item #

600-ALR-AA

NURSE CALL ADAPTER CABLE



Latex-free

Nurse Call Adapter Cable allows the Flex alarm (600-AL2) to plug into a nurse call system.

- 8 feet in length
- Common 1/4" phone plug
- Compatible with most nurse call systems

Item #

600-ALR-AC

PROTECTIVE BOOT



Latex-free

Protective boot protects alarms from accidental drops.

- Durable, long-lasting design
- Attractive teal color

Item #

600-ALPB

HOOK-AND-LOOP STRAP



Hook-and-loop strap gives alarms additional mounting flexibility, easily attaching to wheelchairs and beds.

- Durable, long-lasting design
- 12" in length

Item #

600-ALR-CS

SINGLE CALL CORD



Latex-free

This single call cord is heavy duty and easy to use. It is designed to meet the demanding requirements of institutional healthcare use.

- 8 feet and 10 feet in length
- 1/4" phone plug works with most brands of call systems
- Ergonomically designed pendant is easy to grasp
- Includes security clip to attach pendant end to bed linens

Item #

600-ANCS-8

8 feet in length

600-ANCS-10

10 feet in length

DOUBLE CALL CORD



Latex-free

The double call cord allows two residents to use a single bedside station. Heavy duty and easy to use, this call cord is designed for use in institutional healthcare environments.

- 10 feet in length
- 1/4" phone plug works with most brands of call systems
- Ergonomically designed pendant is easy to grasp
- Includes two security clips to attach each pendant end to bed linens

Item #

600-ANCD-10

GERIATRIC CALL CORD



Latex-free

The geriatric pad is a great alternative for residents who cannot call for assistance using a single call cord. The large pad is easily activated using an entire hand, fist or elbow. Heavy duty and easy to use, this call cord is designed for use in institutional healthcare environments.

- 8 feet in length
- 1/4" phone plug works with most brands of call systems
- Oversized pad is easily activated
- Includes security clip to attach each pendant end to bed linens

Item #

600-ANCS-8G



Professional Medical, Inc.

A Tradition of Quality, Value, & Trusted Service Since 1968

1917 Garnet Ct. New Lenox, IL 60451 ■ 800.648.5190 ■ fax 866.726.7416 ■ promedsupply.com

© Professional Medical, Inc.



We make improving care easier!®

700-Fall Prevention Systems [11/13]