

Pressure Redistribution Surfaces



Proud Member of



PROFESSIONAL MEDICAL, INC.

STATIC PRESSURE REDISTRIBUTION SURFACES

Static pressure redistribution surfaces, including foam mattresses, are non-powered devices that may be indicated when a resident is at risk for pressure ulcer development or delayed healing.

Compliance® Pressure Redistribution Foam Mattress

- Three pressure zones with 670 individual cells
- Antimicrobial, anti-fungal, fluid resistant cover
- Concealed zipper with anti-contamination flap
- WARRANTY: Four years, non-prorated

Item #	Size	Weight Capacity
600-M700-76*	36" x 76" x 6"	300 lbs
600-M700-80*	36" x 80" x 6"	300 lbs
600-M700-84*	36" x 84" x 6"	300 lbs

Embrace® Pressure Redistribution Foam Mattress

- Designated head and foot sections
- Heel section has air channel cut foam for added air circulation and pressure redistribution
- Four-way stretch water-resistant and antimicrobial nylon cover
- Concealed zipper with anti-contamination flap
- WARRANTY: Seven year, non-prorated

Item #	Size	Weight Capacity
600-M800-76*	36" x 76" x 6"	350 lbs
600-M800-80*	36" x 80" x 6"	350 lbs
600-M800-84*	36" x 84" x 6"	350 lbs

Embrace® Elite Pressure Redistribution Gel-Infused Foam Mattress

- Gel infused foam provides unprecedented pressure redistribution, rapid heat dissipation and increased breathability
- Precision cut air pockets allows deeper submersion and increased air circulation
- Lura-Quilt™ aerated foam heel section keeps residents cool and dry
- Heel slope redistributes pressure from the heel to the entire leg
- WARRANTY: 15 year, non-prorated

Item #	Size	Weight Capacity
600-M1000-76*	36" x 76" x 6"	450 lbs
600-M1000-80*	36" x 80" x 6"	450 lbs
600-M1000-84*	36" x 84" x 6"	450 lbs

Embrace® Elite Pressure Redistribution Bariatric Foam Mattress

- Ultra high density, die-cut foam provides exceptional support, comfort and stability
- Lura-Quilt™ aerated foam heel section keeps residents cool and dry
- Four-way stretch water-resistant and antimicrobial nylon cover
- Concealed zipper with anti-contamination flap
- WARRANTY: One year, non-prorated

Item #	Size	Weight Capacity
600-M2000-4276	42" x 76" x 6"	1000 lbs
600-M2000-4280	42" x 80" x 6"	1000 lbs
600-M2000-4284	42" x 84" x 6"	1000 lbs

*Add "RR" to the end of an item number for raised rails

DYNAMIC PRESSURE REDISTRIBUTION SURFACES

Dynamic pressure redistribution surfaces are powered devices that are indicated if a resident can't move into a variety of positions without placing weight on a pressure ulcer, if their weight completely compresses a static device that otherwise maintains its original integrity or if they have a pressure ulcer that is not healing as expected and it is determined that pressure could be contributing to the delay.

Low Air Loss/Alternating Pressure Mattress System

- High-flow pump with 8 L of air per minute
- Adjustable comfort setting allows pressure to be customized according to therapeutic requirements
- Low pressure alarm indicates when the pump's airflow may be compromised and has a mute button
- Static mode provides firm support for egress or ingress
- 18 laser air holes to keep resident cool and dry
- WARRANTY: 18 months

Item #	Size	Weight Capacity
600-AM7000	35" x 80" x 8"	350 lbs

Low Air Loss/Alternating Pressure Mattress System

- High-flow pump with 8 LPM
- Dual therapies low air loss and alternating pressure
- Cell-on-cell air bladder system prevents bottoming out in the event of a power failure
- Transport connector locks in air for easy patient transport
- Quilted top cover provides comfortable sleeping surface
- WARRANTY: 18 months

Item #	Size	Weight Capacity
600-AM8000	35" x 80" x 8"	350 lbs
600-AM8000-RR	35" x 80" x 10"; Raised rails	350 lbs

True Low Air Loss and Alternating Pressure Mattress System

- High-flow pump with 1300 L of air per minute true low air loss
- Dual therapies of low air loss and alternating pressure combines in one pump
- 10 pressure level settings help customize to each individual resident's requirements
- Cell-on-cell air bladder system prevents bottoming out in the event of a power failure
- Auto firm for quick 1 minute mattress inflation
- WARRANTY: 18 months

Item #	Size	Weight Capacity
600-AM9000	36" x 80" x 10"	1000 lbs

Low Air Loss/Alternating Pressure Bariatric Mattress System

- Features alternating pressure cycles of 10, 15, 20 and 25 that allow custom pressure redistribution
- Audible low pressure alarms provide an alert when preset pressure levels decrease
- 10" air cells are constructed of Nylon/PU for added comfort and support
- Pillow top feature prevents unnecessary head and neck movement during alternating mode
- CPR valve provides rapid deflation of mattress
- WARRANTY: Two years

Item #	Size	Weight Capacity
600-AM9000B	42" x 80" x 10"	650 lbs
600-AM9000B48	48" x 80" x 10"	650 lbs





Professional Medical, Inc.

A Tradition of Quality, Value, & Trusted Service Since 1968

We Make Improving Care Easier!®

1917 Garnet Ct. New Lenox, IL 60451

800.648.5190 fax 866.726.7416

promedsupply.com

Compliance®

Improve Care. Save Time. Reduce Costs.



© Professional Medical, Inc.
*We reserve the right to correct any errors
that may occur within this brochure.*

Compliance, Embrace, Smart Choice,
Professional Medical, Inc., ProMed, and
We Make Improving Care Easier! are
trademarks of Professional Medical, Inc.

Orders will be shipped in the order they are received.
Special priced products available while supplies last.
Applicable sales tax will be applied.
Our shipping/handling policy procedures apply.

700-Pressure_Redistribution_Surfaces_Brochure [02/16]