

# Equipment and Furniture Solutions: Spring 2013



PROFESSIONAL MEDICAL, INC.

# Mobility & Lifting



**drive**  
MEDICAL DESIGN & MANUFACTURING



**600-AFLW**

## INNOVATIVE PRODUCT



### LIFTWALKER

- 300 lb weight capacity
- Ideal for individuals who need help sitting or standing
- Easy-to-deploy retractable stand-assist bars
- Simple-to-lock wheels
- WARRANTY: Lifetime limited

Item #	Price
<b>600-AFLW</b>	<b>\$129.95</b>



500-WC16F



500-WC22DF-FA



500-WC18DE



500-WC24DE

SENTRA

SILVER SPORT 2

## STANDARD WHEELCHAIRS

- 250 lb (16"); 300 lb (18"); 350 lb (20") weight capacities
- Powder-coated silver vein steel frame
- Embossed vinyl upholstery
- Durable urethane tires on composite wheels
- WARRANTY: Three-year limited

## OVERSIZED WHEELCHAIRS

- 450 lb weight capacity
- Carbon steel frame with triple-coated chrome
- Double-embossed vinyl upholstery
- Durable urethane tires on composite wheels
- WARRANTY: Three-year limited

## CUSTOMER FAVORITE

### FIXED ARM

Item #	Price
<b>500-WC16F</b>	<b>\$122.95</b>
16" width & swing-away footrests	
500-WC16E	\$132.95
16" width & elevating legrests	
500-WC18F	\$114.95
18" width & swing-away footrests	
500-WC18E	\$122.95
18" width & elevating legrests	

### PADDED DESK ARM - DETACHABLE

Item #	Price
500-WC16DF	\$125.95
16" width & swing-away footrests	
500-WC16DE	\$138.95
16" width & elevating legrests	
500-WC18DF	\$125.95
18" width & swing-away footrests	
<b>500-WC18DE</b>	<b>\$138.95</b>
18" width & elevating legrests	
500-WC20DF	\$144.95
20" width & swing-away footrests	
500-WC20DE	\$149.95
20" width & elevating legrests	
500-WC22DF	\$264.10
22" width & swing-away footrests	
500-WC22DE	\$269.25
22" width & elevating legrests	
500-WC24DF	\$319.95
24" width & swing-away footrests	
<b>500-WC24DE</b>	<b>\$329.95</b>
24" width & elevating legrests	

### FULL ARM - DETACHABLE

Item #	Price	Item #	Price
500-WC16DF-FA	\$125.95	500-WC20DE-FA	\$159.95
16" width & swing-away footrests		20" width & elevating legrests	
500-WC16DE-FA	\$138.95	<b>500-WC22DF-FA</b>	<b>\$264.10</b>
16" width & elevating legrests		22" width & swing-away footrests	
500-WC18DF-FA	\$138.95	500-WC22DE-FA	\$269.25
18" width & swing-away footrests		22" width & elevating legrests	
500-WC18DE-FA	\$132.95	500-WC24DF-FA	\$319.95
18" width & elevating legrests		24" width & swing-away footrests	
500-WC20DF-FA	\$154.95	500-WC24DE-FA	\$329.95
20" width & swing-away footrests		24" width & elevating legrests	

**CHECK OUT**  
promedsupply.com

for full descriptions  
on these and our  
other items



500-GC-3PM



500-WC18DE-REC



500-WC30DE  
FOOT AND LEGRESTS  
SOLD SEPARATELY

## SENTRA



500-STDS860



500-STDS803



500-STDS802



600-BG60-SS

600-BG60-P



### 3-POSITION GERI RECLINER CHAIR

- 250 lb weight capacity
- Three reclining positions
- **WARRANTY:** *Frame:* Lifetime limited; *Upholstery and moving parts:* One year

Item #	Price
<b>500-GC-3PM</b> Mauve	\$365.95
500-GC-3PBR Blue	\$365.95
500-GC-3PJ Jade	\$365.95

### RECLINING WHEELCHAIRS

- 300 lb (16"/18"); 350 lb (20"); 450 lb (22") weight capacities
- State-of-the-art hydraulic reclining mechanism for 180° adjustments
- Composite mag-style wheels are light and maintenance-free
- Padded armrests for resident comfort
- Wheels set back on frame to prevent tipping
- **WARRANTY:** Lifetime limited

Item #	Price
500-WC16DE-REC 16" detachable desk arm; elevating legrests	\$298.95
<b>500-WC18DE-REC</b> 18" detachable desk arm; elevating legrests	\$298.95
500-WC18FE-REC 18" detachable full arm; elevating legrests	\$298.95
500-WC20DE-REC 20" detachable desk arm; elevating legrests	\$298.95
500-WC22DE-REC 22" detachable desk arm; elevating legrests	\$418.95

### BARIATRIC WHEELCHAIRS

- 700 lb weight capacity
  - Carbon steel frame with triple-coated chrome for a chip-proof finish
  - Reinforced steel gussets at all weight-bearing points
  - Extra-large padded armrests
  - **WARRANTY:** Lifetime limited
- Select footrest (500-PHSF) or elevating legrests (500-PHEL) when ordering*

Item #	Price
500-WC26DE 26" detachable desk arm	\$598.95
500-WC28DE 28" detachable desk arm	\$655.80
<b>500-WC30DE</b> 30" detachable desk arm	\$689.95
500-WC26DE-FA 26" detachable full arm	\$761.54
500-WC28DE-FA 28" detachable full arm	\$655.80
500-WC30DE-FA 30" detachable full arm	\$666.80
500-PHSF 26"-30" swingaway footrest	\$76.50 (PR)
500-PHEL 26"-30" elevating legrests	\$82.50 (PR)

### REPLACEMENT PARTS AND ACCESSORIES

Item #	Price
<b>500-STDS860</b> Front anti-tippers	\$55.05 (PR)
<b>500-STDS802</b> Back anti-tippers with wheels	\$38.74 (PR)
<b>500-STDS803</b> Oxygen "e" cylinder holder	\$35.90 (EA)
500-WCARMS-D Replacement desk length black armrests	\$11.86 (EA)
500-WCARMS-F Replacement full length black armrests	\$11.86 (EA)
<b>600-BG60-P</b> Gait belt with plastic buckle	\$7.35
<b>600-BG60-SS</b> Gait belt with stars and stripes	\$9.45

### FOOTRESTS AND ELEVATING LEGRESTS

Item #	Price
500-STDDSF 16"-20" swingaway footrest	\$35.74 (PR)
500-STDEL 16"-20" elevating legrests	\$48.45 (PR)
500-STDSF 22"-24" swingaway footrest	\$35.95 (PR)
500-STDEL 22"-24" elevating legrests	\$48.45 (PR)



**500-LP350-S**



**500-LP600\***



**500-LP450-L\***



**500-S45L-S**



**500-S355**

\*slings not included

## SCHEDULE A DEMO

Contact your territory  
manager or call us at  
**(800) 648-5190**

### STAND-ASSIST LIFT

- 400 lb weight capacity
- Boom height: 26.2" to 81.1"
- Base width: 24.5" closed; 40.2" open
- 4.5" base clearance
- Easy-roll 4" casters
- Four sling attachment points for a variety of stand-up sling designs
- **WARRANTY:** Frame: Lifetime; Electronic components and actuator: One year

Item #	Price
<b>500-LP350-S</b>	<b>\$2,495.99</b>
<b>500-S355</b>	<b>\$99.99</b>
Standing sling; 440 lb capacity	

### BARIATRIC BATTERY- POWERED LIFT

- 600 lb weight capacity
- Boom height: 32" to 74.5"
- Base width: 25.75" closed; 45.75" open
- 2" base clearance
- Six-point cradle
- Extra-wide coated cradle with 360° rotation
- **WARRANTY:** Frame: Lifetime; Electronic components and actuator: One year

Item #	Price
<b>500-LP600</b>	<b>\$2,039.95</b>
600 lb capacity	

### BATTERY-POWERED LIFT

- 450 lb weight capacity
- Boom height: 27.5" to 76"
- Base width: 15" closed; 41" open
- Low 2" base clearance
- For use with two-point and four-point slings
- Fast, easy charge from AC outlet
- **WARRANTY:** Lifetime limited

Item #	Price
<b>500-LP450</b>	<b>\$1,385.99</b>
450 lb capacity	
<b>500-LP450-L</b>	<b>\$1,579.48</b>
With removable, rechargeable battery; 450 lb capacity	

### FOUR POINT AND STANDING SLINGS


Item #	Price
<b>500-S45L-S</b>	<b>\$105.29</b>
600 lb capacity; Large	
<b>500-S45M-S</b>	<b>\$105.29</b>
600 lb capacity; Medium	
<b>500-S45L-MC</b>	<b>\$89.50</b>
Mesh with commode; 600 lb capacity; Large	
<b>500-S45M-MC</b>	<b>\$89.50</b>
Mesh with commode; 600 lb capacity; Medium	

## CUSTOMER FAVORITE

### APEX GEL-FOAM CUSHION

- High-performance fabric promotes healing via moisture evaporation and cooling at pressure points
- Dual chamber gel-pod reduces gel migration and eliminates bottoming out
- Breathable four-way stretch removable cover is antimicrobial and mildew-resistant with an anti-slip bottom
- Easy-to-attach quick-release straps
- Meets CA 117 standards

Item #	Price
<b>600-WC1816DGF-EA</b> 18" x 16" x 2"	<b>\$59.95</b>
600-WC2016DGF-EA 20" x 16" x 2"	<b>\$79.95</b>



**600-WC1816DGF-EA**

We offer a comprehensive line of cushions and accessories designed to help protect against skin breakdown, improve resident comfort, correct positioning and reduce sliding. The innovative APEX designs feature dual-density foam layers for pressure redistribution and a dual gel-pod for increased comfort and stability. These high-quality items are proudly made in the USA.



**SHOP NOW**  
promedsupply.com

for these and other support surfaces

## CUSTOMER FAVORITE

### APEX WEDGE GEL-FOAM CUSHION



- Wedge shape for improved posture and slide control
- Breathable four-way stretch cover promotes active healing
- Slip-resistant bottom with quick-release straps
- Meets CA 117 standards

Item #	Price
<b>600-WC1816WGF-EA</b> 18" x 16" x 4" to 2"	<b>\$69.95</b>
<b>600-WC2016WGF-EA</b> 20" x 16" x 4" to 2"	<b>\$89.95</b>


### APEX LEG ABDUCTOR (POMMEL) GEL-FOAM CUSHION



- Contoured foam pommel for comfortable abduction
- Firm foundation eliminates hammocking
- Slip-resistant bottom with quick-release straps
- Meets CA 117 standards

Item #	Price
<b>600-WC1816PEA</b> 18" x 16" x 3"	<b>\$79.50</b>

### GEL-FOAM CUSHION



- High-density foam and a dual-chamber gel pod provide superior pressure relief and comfort
- Low shear, antimicrobial, fluid-proof, fire-resistant
- Slip-resistant bottom with Velcro straps to keep cushion in place
- Meets CA 117 standards

Item #	Price
<b>600-WC1816D-EA</b> 18" x 16" x 2"	<b>\$22.95</b>
<b>600-WC1816-EA</b> 18" x 16" x 3"	<b>\$26.95</b>
<b>600-WC2016-EA</b> 20" x 16" x 3"	<b>\$39.50</b>
<b>600-WC2218-EA</b> 22" x 18" x 3"	<b>\$44.50</b>
<b>600-WC2418-EA</b> 24" x 18" x 3"	<b>\$52.50</b>

### WHEELCHAIR CUSHION



- Combines a foam base and chambered gel insert for enhanced comfort and pressure relief
- Low-shear, antimicrobial, fire-resistant and fluid-proof removable cover
- Slip-resistant bottom with Velcro straps to keep cushion in place
- Meets CA 117 standards

Item #	Price
<b>500-WC02V</b> 16" x 18" x 2"	<b>\$16.95</b>
<b>500-WC03V</b> 16" x 18" x 3"	<b>\$16.95</b>
<b>500-WC04V</b> 16" x 18" x 4"	<b>\$17.95</b>

### HEEL SUSPENSION BOOT



- Super-soft antimicrobial peachskin cloth protects thin, fragile skin
- Hook and loop closures allow boot to be fastened in multiple positions
- Extended, cushioned design distributes heel weight over a larger area, reducing pressure and increasing blood flow

Item #	Price
600-PHE3-M Medium	\$21.95
600-PHE3-L Large/Extra Large	\$26.95

### HEEL FLOTATION BOOT



- Soft velour fabric
- Adaptable hook and loop closures
- Water-based gel pack for hot or cold therapy
- One size fits most (*average adult*)

Item #	Price
600-PHG3	\$23.95

## LAP POSITIONERS



600-WCLH

600-WCLH-D

600-WCLH-F

### LAP POSITIONERS

- Recommended for residents needing upper torso support or an elevated resting spot for hands
- Assist with upper-body alignment and increase safety by gently reminding residents to stay seated
- Available in self releasing, full arm and desk arm styles

Item #	Price
<b>600-WCLH</b> Self releasing; 18-22" x 20" x 4"	\$31.95
<b>600-WCLH-F</b> Full arm; 27" x 18" x 3"	\$31.95
<b>600-WCLH-D</b> Desk arm; 27" x 18" x 3"	\$31.95

### BALANCED AIRE CUSHION™



- 300 lb weight capacity
- Made from flame-resistant neoprene
- Includes a manual inflation pump, patch kit and breathable, low-shear cover
- **WARRANTY:** *Cushion material and workmanship:* 18 months; *Exterior cover:* 30 days

Item #	Price
500-WC1816-2LP 18"x16"x 2"	\$245.95
500-WC1816-4HP 18"x16"x 4"	\$245.95

### 2-BUTTON FOLDING WALKER AND ACCESSORIES



- 350 lb weight capacity
- Adjustable height: 32"-39"
- Easy push-button mechanisms
- U-shaped frame for greater clearance
- **WARRANTY:** Lifetime limited

Item #	Price
500-W101	\$27.95
500-W101-5WF 5" fixed wheels ( <i>shown above</i> )	\$30.25
500-W1A-5WF 5" fixed wheels with rear cap glides	\$10.78 (PR)
500-W1A-GC Cap glides	\$5.30 (PR)
500-W1A-GT Tennis ball glides	\$14.92 (PR)

### FOLDING STEEL COMMODE



- 350 lb weight capacity
- Folds flat for convenient storage and transportation
- Durable plastic snap-on seat and lid
- Includes a 7.5-quart commode bucket, cover and splash shield
- **WARRANTY:** Lifetime limited

Item #	Price
500-C101	\$48.19
500-C101-P 12-qt replacement pail and cover	\$5.50

### BIARIATRIC FOLDING COMMODE



- 650 lb weight capacity
- Folds flat for convenient storage and transportation
- Removable tool-free back
- Durable 16.5" snap-on seat with lid
- Includes a 12-quart commode bucket with handle, cover and splash shield
- **WARRANTY:** Lifetime limited

Item #	Price
500-C101-XW	\$99.99
500-C101-P 12-qt replacement pail and cover	\$5.50

# Beds & Accessories



**CHECK OUT**  
[promedsupply.com](http://promedsupply.com)

for full descriptions on  
these and other items

#### BASIC HIGH-LOW BED

- 450 lb weight capacity
- 80" L x 35" W
- Height range: 8" to 26"
- Roll-at-any-height design
- Locking system
- WARRANTY: *Frame:* 15 years;  
*Electronic components:* Two years;  
*Welds:* Lifetime

Item #	Price
<b>600-BED1</b>	<b>\$998.65</b>



**600-BED1**



## HIGH-LOW BEDS



600-BED2



600-BED3

Smokey Hollow  
DEFAULT SHAPE

SOLAR  
OAK

AMBER  
CHERRY

Mill Creek

SOLAR  
OAK

AMBER  
CHERRY

### SMOKEY HOLLOW HEADBOARDS AND FOOTBOARDS

Item #	Price
600-BED-HB1 Headboard; Rich Solar Oak laminate	\$42.60
600-BED-FB1 Footboard; Rich Solar Oak laminate	\$42.60
600-BED-HB2 Headboard; Rich Amber Cherry laminate	\$42.60
600-BED-FB2 Footboard; Rich Amber Cherry laminate	\$42.60

### MILL CREEK HEADBOARDS AND FOOTBOARDS

Item #	Price
600-BED-HB3 Headboard; Rich Solar Oak laminate	\$42.60
600-BED-FB3 Footboard; Rich Solar Oak laminate	\$42.60
600-BED-HB4 Headboard; Rich Amber Cherry laminate	\$42.60
600-BED-FB4 Footboard; Rich Amber Cherry laminate	\$42.60

### BASIC+ HIGH-LOW BED

- 450 lb weight capacity
- 80" L x 35" W
- Height range: 8" to 26"
- Advanced caster locking system
- Optional 42" expandable deck
- WARRANTY: *Frame:* 15 years;  
*Electronic components:* Two years;  
*Welds:* Lifetime

Item #	Price
600-BED2	\$1,138.65

### ULTIMATE HIGH-LOW BED




- 500 lb weight capacity
- Expands from 76" to 84" L and 35" to 42" W
- Height range: 6.88" to 30.88"
- Roll-at-any-height design
- Provides more workable surface for repositioning and care
- WARRANTY: *Frame:* 15 years;  
*Electronic components:* Two years;  
*Welds:* Lifetime

Item #	Price
600-BED3	\$1,618.65

# Beds & Accessories



## AVAILABLE COLORS

- WILD CHERRY 
- SOLAR OAK 
- KENSINGTON MAPLE 

## LIBERTY BED

- 450 lb weight capacity
- 76" or 80" L x 35" W
- Adjustable brackets allow bed to be used at both 76" and 80" lengths
- Height range: 9.9" to 24.4"
- 360-degree mobility at any height
- Quiet, smooth operation with three DC motors
- Surface treated with antimicrobial SilverSolutions
- WARRANTY: Three years

Item #	Price
075-LBT7174	\$735.95
Liberty Bed (assist rails sold separately)	
075-AB53600	\$59.95
Headboard & footboard; Solar Oak	
500-BSC3-O	\$415.99
3-drawer bedside cabinet; Solar Oak	
500-DRS3-O	\$523.99
4-drawer chest; Solar Oak	

## AUTO-TOUCH OVERBED TABLE



- Adjustable from 28" to 45"
- 15" x 30" tabletop
- Walnut wood grain laminated top
- 2" swivel casters for easy maneuverability
- WARRANTY: Lifetime limited

Item #	Price
500-TTBED	\$49.99

## SAFETY HALF-LENGTH RAILS



- Fits all spring-style healthcare beds
- No-gap style
- Steel construction
- WARRANTY: Lifetime limited

Item #	Price
500-RAILH	\$56.95

## HALF AND FULL VINYL BED RAIL PADS



- 34" or 69" L x 11" W x 1" D
- Half-length and full-length styles available
- Provide protective padding to reduce the risk of resident injury

Item #	Price
600-BRP1134-PR	\$21.25
Half-length vinyl bed rail pads; rails sold separately	
600-BRP1169-PR	\$37.95
Full-length vinyl bed rail pads; rails sold separately	

## BARIATRIC TRAPEZE BAR



- 750 lb weight capacity
- Heavy-duty grab handle and multi-height position strap
- Easy-to-use rung-style strap attachments for quick relocation
- Freestanding design




Item #	Price
600-T01B-F	\$395.55
Freestanding	

**ZENITH 9000  
BED**

- 600 lb weight capacity
- 76" or 80" L x 35" W
- Adjustable brackets allow bed to be used at both 76" and 80" lengths
- Height range: 7" to 30"
- 360-degree mobility at any height
- Fast, even-speed Timesaver Technology™ motors
- Surface treated with antimicrobial SilverSolutions
- **WARRANTY:** *Frame:* 15 years; *Motors/control box:* Five years; *Other electronics:* Three years

Item #	Price
075-ZZ98174	\$1,279.95
Assist rails sold separately	
075-Z04500CCR	\$61.25
Headboard & footboard; Wild Cherry	
500-BSC3-C	\$415.99
3-drawer bedside cabinet; Wild Cherry	
500-DRS3-C	\$523.99
4-drawer chest; Wild Cherry	

**AVAILABLE COLORS**

- WILD CHERRY** 
- SOLAR OAK** 
- KENSINGTON MAPLE** 



**BARIATRIC**




**FREE  
42" FOAM MATTRESS  
AND PAIR OF "T" RAILS**

**CUSTOMER FAVORITE**

**FULL ELECTRIC  
BARIATRIC BED**

- 600 lb weight capacity
- 88" L x 42" W
- Height range: 18.25" to 26.25"
- Heavy-duty frame for strength and resident safety
- Split plan design with removable bed ends for easy setup
- Quiet, smooth operation
- **WARRANTY:** *Frame:* Five years; *Other parts and components:* One year; *Welds:* Lifetime

Item #	Price
500-BB600-PKG Package	\$2,249.00



**Smokey Hollow**  
DEFAULT SHAPE

**SOLAR OAK**

**BARIATRIC  
EXPANDABLE BED**

- 750 lb weight capacity
- 80"-88" L x 39"-48" W
- Height range: 15.75" to 29.25"
- Smokey Hollow headboard and footboard available in Rich Solar Oak laminate
- **WARRANTY:** *Frame:* 15 years; *Electronic components:* Two years; *Welds:* Lifetime

Item #	Price
600-BED4	\$3,149.65
600-BED4-HB1 Headboard; Rich Solar Oak laminate	\$77.50
600-BED4-FB1 Footboard; Rich Solar Oak laminate	\$77.50



500-BSC1-C






500-BSC2-C

500-DRS2-C



AVAILABLE COLORS  
HIGH-QUALITY LAMINATE

- CHERRY 
- OKLAHOMA OAK 
- WHITE ROCK MAPLE 



500-WRD1-2-C

**SHOP NOW**  
promedsupply.com

for these and other furniture options

BESIDE CABINETS

- 29" W x 18" D x 29.5" H
- Durable nickel-finish hardware
- All rounded corners for safety
- WARRANTY: Five-year limited

Item #	Price
<b>500-BSC1-C</b>	<b>\$129.99</b>
1-door, 1-drawer cabinet; Cherry	
500-BSC1-O	\$129.99
1-door, 1-drawer cabinet; Oklahoma Oak	
500-BSC1-M	\$129.99
1-door, 1-drawer cabinet; White Rock Maple	
<b>500-BSC2-C</b>	<b>\$138.99</b>
3-drawer cabinet; Cherry	
500-BSC2-O	\$138.99
3-drawer cabinet; Oklahoma Oak	
500-BSC2-M	\$138.99
3-drawer cabinet; White Rock Maple	

3-DRAWER AND 4-DRAWER CHESTS

- 29" W x 18" D x 29.5" H (3-dr.); 29" W x 18" D x 37" H (4-dr.)
- Durable nickel-finish hardware
- All rounded corners for safety
- WARRANTY: Five-year limited

Item #	Price
500-DRS1-C	\$211.99
3-drawer chest; Cherry	
500-DRS1-O	\$211.99
3-drawer chest; Oklahoma Oak	
500-DRS1-M	\$211.99
3-drawer chest; Maple	
<b>500-DRS2-C</b>	<b>\$222.99</b>
4-drawer chest; Cherry	
500-DRS2-O	\$222.99
4-drawer chest; Oklahoma Oak	
500-DRS2-M	\$222.99
4-drawer chest; White Rock Maple	

1-DOOR, 2-DRAWER WARDROBE

- 22" W x 23" D x 78" H
- Durable nickel-finish hardware
- All rounded corners for safety
- WARRANTY: Five-year limited

Item #	Price
<b>500-WRD1-2-C</b>	<b>\$386.99</b>
1-door, 2-drawer wardrobe; Cherry	
500-WRD1-2-O	\$386.99
1-door, 2-drawer wardrobe; Oklahoma Oak	
500-WRD1-2-M	\$386.99
1-door, 2-drawer wardrobe; White Rock Maple	

**500-WRD2-C**



**500-WRD2-2-C**



**2-DOOR  
WARDROBE**

- 30" W x 23" D x 78" H
- Durable nickel-finish hardware
- All rounded corners for safety
- WARRANTY: Five-year limited

Item #	Price
<b>500-WRD2-C</b> 2-door wardrobe; Cherry	<b>\$442.99</b>
500-WRD2-O 2-door wardrobe; Oklahoma Oak	<b>\$442.99</b>
500-WRD2-M 2-door wardrobe; White Rock Maple	<b>\$442.99</b>

**2-DOOR, 2-DRAWER  
WARDROBE**

- 30" W x 23" D x 78" H
- Durable nickel-finish hardware
- All rounded corners for safety
- WARRANTY: Five-year limited

Item #	Price
<b>500-WRD2-2-C</b> 2-door, 2-drawer wardrobe; Cherry	<b>\$484.99</b>
500-WRD2-2-O 2-door, 2-drawer wardrobe; Oklahoma Oak	<b>\$484.99</b>
500-WRD2-2-M 2-door, 2-drawer wardrobe; White Rock Maple	<b>\$484.99</b>



The **Victory Blood Glucose Meter**, our Compliance Solution, is a technologically advanced auto-code meter that delivers accurate results in just five seconds. The accompanying test strips feature gold-plated electrodes for exceptionally consistent readings.

**Features:**

- Test strip release button protects users from direct contact with blood
- Small 0.5 µl sample size
- Test range: 20-600mg/dl
- Glucose dehydrogenase enables the use of capillary, venous and arterial blood for glucose measurement
- Allows for alternate site testing

ITEM#	DESCRIPTION
600-CBG1	Blood glucose meter (ea)
600-CBG2	High-low control solutions (st)
600-CBG50	Blood glucose test strips (50/bx 10bx/ctn)
600-P200A	Alcohol prep pads (200/bx 20bx/cs)



**VICTORY**  
**Blood Glucose System**

Designed with long-term care in mind



600-M700-76



**FREE SHIPPING  
WHEN YOU BUY 5**

Qualify by buying 5 or more foam mattresses of any size or any combination


### 3-YEAR PRESSURE REDISTRIBUTION MATTRESS

- 300 lb weight capacity
- Three pressure zones
- Channel-cut high-density polyurethane foam
- **WARRANTY:** Three years, non-prorated

Item #	Price
<b>600-M700-76</b> 36" x 76" x 6"	<b>\$139.95</b>
600-M700-80 36" x 80" x 6"	<b>\$139.95</b>
600-M700-84 36" x 84" x 6"	<b>\$149.95</b>

## CUSTOMER FAVORITE

### 6-YEAR PRESSURE REDISTRIBUTION MATTRESS



- 350 lb weight capacity
- 4-way stretch nylon waterproof cover
- Provides additional pressure reduction in heel area
- Solid-edge support for easy resident transfer
- **WARRANTY:** Six years

Item #	Price
600-M850-76 36" x 76" x 6"	<b>\$187.95</b>
600-M850-80 36" x 80" x 6"	<b>\$187.95</b>
600-M850-84 36" x 84" x 6"	<b>\$198.95</b>

### 15-YEAR PRESSURE REDISTRIBUTION MATTRESS



- 450 lb weight capacity
- 4-way stretch nylon waterproof cover
- Designated head and foot sections
- Concealed zipper with Barrier Stop Over-Flap
- **WARRANTY:** 15 years

Item #	Price
600-M900-76 36" x 76" x 6"	<b>\$214.99</b>
600-M900-80 36" x 80" x 6"	<b>\$214.99</b>
600-M900-84 36" x 84" x 6"	<b>\$225.99</b>


### DYNAMIC ELITE MATTRESS



- 400 lb weight capacity
- Quickly channels heat away from the body, helping to reduce pressure ulcer incidence
- Provides superior pressure redistribution
- **WARRANTY:** 10 years, non-prorated

Item #	Price
600-M1000-76 36" x 76" x 6"	<b>\$229.99</b>
600-M1000-80 36" x 80" x 6"	<b>\$229.99</b>
600-M1000-84 36" x 84" x 6"	<b>\$245.99</b>

### DYNAMIC ELITE MATTRESS WITH RAISED RAILS



- 400 lb weight capacity
- Quickly channels heat away from the body, helping to reduce pressure ulcer incidence
- Provides superior pressure redistribution
- Side rails decrease fall risk
- **WARRANTY:** 10 years, non-prorated

Item #	Price
600-M1000-76RR 36" x 76" x 6"	<b>\$255.99</b>
600-M1000-80RR 36" x 80" x 6"	<b>\$264.99</b>
600-M1000-84RR 36" x 84" x 6"	<b>\$269.99</b>

## CUSTOMER FAVORITE

### ALTERNATING PRESSURE MATTRESS SYSTEM



- 350 lb weight capacity
- Includes 20-bladder mattress and high-flow pump
- Helps prevent and treat pressure sores
- Adjustable comfort settings allow pressure to be customized by resident weight
- WARRANTY: Pump: One year

Item #	Price
600-AM7000 36" x 80" x 8"	\$465.95

### LOW AIR LOSS & ALTERNATING PRESSURE MATTRESS SYSTEM



- 350 lb weight capacity
- Provides low air loss and alternating pressure therapy in 10, 15, 20 and 25-minute cycles
- 10 adjustable comfort cycles allow pressure to be customized by resident weight
- Pillow top keeps head and neck stable
- WARRANTY: Limited

Item #	Price
600-AM8000 36" x 80" x 8"	\$789.95

### TRUE LOW AIR LOSS MATTRESS SYSTEM WITH PULSATE



- 1000 lb weight capacity
- Includes mattress and LAL/pulsate pump
- Provides a cool and dry sleeping surface
- Exclusive pulsate mode reduces pressure and increases blood flow
- WARRANTY: Limited

Item #	Price
600-AM9000 36" x 80" x 10"	\$2,015.99

### BARIATRIC LOW AIR LOSS WITH ALTERNATING PRESSURE MATTRESS



- 1000 lb weight capacity
- Provides low air loss and alternating pressure therapies
- Microprocessor unit provides continuous feedback to equalize pressure and maximize comfort
- WARRANTY: Limited

Item #	Price
600-AM9000B 42" x 80" x 10"	\$2,189.99
600-AM9000B48 48" x 80" x 10"	\$2,209.99



## New Equipment & Furniture Now

*Tailored leasing plans to fit your needs and budget*

Professional Medical's exclusive **Loan-To-Own (LTO) Program** can bolster your facility's therapeutic equipment inventory, help you acquire new state-of-the-art merchandise and save you thousands of dollars over other traditional rental programs.

**Use your funds wisely now, not during survey time!** Allow Professional Medical to help evaluate the equipment needed to provide the safety your employees expect and the quality care your residents deserve. Now, with our LTO solution, you can afford the equipment needed to provide top-notch care while remaining compliant with state and federal guidelines.

Your ProMed sales representative can walk you through preparing a customized budget analysis offering the best therapeutic equipment backed by the most cost-effective leasing program in the industry.





# DESIGNING BETTER CARE





## When it's time to renovate your facility, turn to the experts.

Professional Medical has partnered with Basic American to create a complimentary design program that aims to enhance resident comfort and wellness. We want to help you create a warm, home-like environment that your residents and their families will enjoy.

Not only will you avoid having to pay an interior designer's hourly charges, you'll also be partnering with a team of professionals who are in tune with the specific challenges of designing for long-term care. Our team is detail-oriented and aware of the codes and regulations that need to be taken into account when designing a space that is both attractive and functional.

We provide the following services, all at no charge:

- Space planning
- Schematic design
- Design development
- Furniture selection
- Color and materials selection
- Custom furniture and millwork design
- Artwork selection
- Graphics and signage
- Accessories and lighting specifications
- Installation administration
- In-house purchasing




**AFTER**



**Before**

Whether you're looking to spruce up or completely overhaul resident rooms, dining areas, activity spaces, common areas or office settings, our experienced interior designers are ready to help you create the perfect look for your facility.

**To learn more, contact your ProMed territory manager or call our Customer Care team at (800) 648-5190.**



# Shower Chairs & Hydration Carts



500-PVC-SC18




500-PVC-SC18FF



500-PVC-SC26

## STANDARD MESH AND SOLID VINYL COLORS - SPECIFY UPON ORDER



Our shower chairs and hydration/distribution carts are made from  **durable, healthcare-grade PVC**  that inhibits the growth of bacteria and is free of phthalates and lead.

These products are reinforced at all stress-related areas for increased durability and performance. They won't rust, fade or stain!

### PVC SHOWER CHAIR WITH CASTERS

- 300 lb weight capacity
- Internal width: 18"; external width: 22"
- 3" or 5" casters
- Push-pull handle for easy maneuvering
- Fast-drying removable mesh sling
- WARRANTY: PVC frame: Five-year limited

Item #	Price
<i>Specify mesh color</i>	
<b>500-PVC-SC18</b>	<b>\$99.95</b>
3" casters	
500-PVC-SC18-C5	<b>\$136.95</b>
5" casters	

### PVC SHOWER CHAIR WITH FOLDING FOOTREST

- 300 lb weight capacity
- Internal width: 18"; external width: 22"
- 3" twin rust-proof casters
- Folding footrest
- Fast-drying removable mesh sling
- WARRANTY: PVC frame: Five-year limited

Item #	Price
<i>Specify mesh color</i>	
<b>500-PVC-SC18FF</b>	<b>\$135.95</b>

### PVC BARIATRIC SHOWER CHAIR

- 425 lb weight capacity
- Internal width: 26"; external width: 30"
- 4" rust-proof casters
- Push-pull handle for easy maneuvering
- Fast-drying removable mesh sling
- WARRANTY: PVC frame: Five-year limited

Item #	Price
<i>Specify mesh color</i>	
<b>500-PVC-SC26</b>	<b>\$195.95</b>



**500-PVC-HC**



**500-PVC-HC2**

**36-QUART  
HYDRATION CART**

- 36-quart ice chest
- Ice scoop
- Storage shelf
- 3" twin rust-proof casters
- WARRANTY: PVC frame: Five-year limited

Item #	Price
<i>Specify mesh or solid vinyl color</i>	
<b>500-PVC-HC</b>	<b>\$192.99</b>

**48-QUART  
HYDRATION CART**

- 48-quart ice chest
- Ice scoop
- Storage shelf
- Pull-out storage/drain tub
- 3" twin rust-proof casters
- WARRANTY: PVC frame: Five-year limited

Item #	Price
<i>Specify mesh or solid vinyl color</i>	
<b>500-PVC-HC2</b>	<b>\$214.99</b>

Our **fixed height bath** is ideal for long-term care use, especially when transfer from a wheelchair is required. The powered traverse makes it easy for residents to get in and out of the bath, helping them retain their independence. The seat can be easily lowered outside the tub to a convenient height for comfort and safety.

The bath also features high-end touches such as a stylish, extra-deep tub and a convenient hand control. Battery back-up comes standard for total peace of mind.

**Additional features:**

- 330 lb lifting capacity
- 60" and 66" lengths available
- Left- and right-hand models available
- Seat can be set at any height
- Smooth, easy-clean surface
- WARRANTY: Two years

Item #	Price
500-TUB1	<b>\$10,497.65</b>



# THE FIXED HEIGHT BATH

# Hampers & Laundry Carts



500-PVC-HM55-1

500-PVC-HM55-2



500-PVC-LC25-3

500-PVC-LC25-3A

500-PVC-LC67-6A



**Innovative Product**

## HAMPER WITH MESH ZIPPER BAG

- Zipper-front opening eliminates lifting heavy loads over the top of the hamper
- Durable high-impact plastic lid
- 3" deluxe, heavy-duty reinforced casters
- Floor plate to support retainer bag
- WARRANTY: PVC frame: Five-year limited

Item #	Price
<i>Specify mesh or solid vinyl color</i>	
<b>500-PVC-HM55-1</b>	<b>\$129.99</b>
Single hamper; 26 gallon	
<b>500-PVC-HM55-2</b>	<b>\$210.99</b>
Double hamper; 52-gallon	

## SOLID SHELF LINEN CART

- 100 lb carrying capacity per shelf
- Shelves are 12" H x 20" W (see below for available lengths)
- 10" W handles
- WARRANTY: PVC frame: Five-year limited

Item #	Price
<i>Specify mesh or solid vinyl color</i>	
500-PVC-LC25-2	\$154.99
2-shelf cart; 20" x 25"	
500-PVC-LC32-2	\$194.86
2-shelf cart; 20" x 32"	
500-PVC-LC45-2	\$252.99
2-shelf cart; 20" x 45"	
500-PVC-LC50-2	\$268.99
2-shelf cart; 20" x 50"	
<b>500-PVC-LC25-3</b>	<b>\$183.72</b>
3-shelf cart; 20" x 25"	
500-PVC-LC32-3	\$254.99
3-shelf cart; 20" x 32"	
500-PVC-LC45-3	\$278.29
3-shelf cart; 20" x 45"	
500-PVC-LC50-3	\$342.99
3-shelf cart; 20" x 50"	

## PERFORATED SHELF LINEN CART

- 75 lb carrying capacity per shelf
- Shelves are 12" H x 20" W (see below for available lengths)
- 10" W handles
- WARRANTY: PVC frame: Five-year limited

Item #	Price
<i>Specify mesh or solid vinyl color</i>	
500-PVC-LC25-2A	\$154.99
2-shelf cart; 20" x 25"	
500-PVC-LC32-2A	\$194.86
2-shelf cart; 20" x 32"	
500-PVC-LC45-2A	\$243.58
2-shelf cart; 20" x 45"	
500-PVC-LC50-2A	\$277.84
2-shelf cart; 20" x 50"	
<b>500-PVC-LC25-3A</b>	<b>\$200.38</b>
3-shelf cart; 20" x 25"	
500-PVC-LC32-3A	\$242.04
3-shelf cart; 20" x 32"	
500-PVC-LC45-3A	\$340.21
3-shelf cart; 20" x 45"	
500-PVC-LC50-3A	\$342.27
3-shelf cart; 20" x 50"	

## DISTRIBUTION CART

- Shelf: 20" W x 43.63" L
- Length with handles: 50"; height with casters: 66.38"
- Hanging rack: 36.75"
- Bin: 16" W x 12.25" L x 5.5" D
- WARRANTY: PVC frame: Five-year limited

Item #	Price
<i>Specify mesh or solid vinyl color</i>	
<b>500-PVC-LC67-6A</b>	<b>\$338.99</b>



**500-HM-1**



**500-HM-2**



**500-HM-3**



**500-3PPS-A**

**SINGLE "EASY ACCESS" HAMPER WITH FOOT PEDAL**

- 7/8" high-quality tubular steel with a chrome finish
- Includes one Easy Access bag
- Foot pedal
- 4" non-marking, quiet-ride casters (two locking)
- WARRANTY: Limited

Item #	Price
<b>500-HM-1</b>	<b>\$219.99</b>

**SINGLE MEDIUM-DUTY HAMPER WITH FOOT PEDAL**

- 7/8" high-quality tubular steel with a chrome finish
- Bag lock secures liners in place
- Foot pedal
- 3" non-marking, quiet-ride casters (two locking)
- WARRANTY: Limited

Item #	Price
<b>500-HM-2</b>	<b>\$81.99</b>

**3-BUSHEL POLY LAUNDRY HAMPER**

- Poly construction with sturdy zinc-plated base
- Leakproof and virtually indestructible
- Foot pedal
- 2" hooded ball casters for smooth transporting
- WARRANTY: Limited

Item #	Price
<b>500-HM-3</b>	<b>\$157.99</b>

**HAMPER LIDS**

Item #	Price
500-HML White	\$26.99
500-HMLID-BLU Blue	\$26.99
500-HMLID-BLK Black	\$26.99
500-HMLID-RED Red	\$26.99
500-HMLID-YEL Yellow	\$26.99

**HAMPER LABELS**

Item #	Price
500-HMLAB-BLU Blue, "Soiled Linen" lettering	\$36.99
500-HMLAB-WHT White, "Soiled Linen" lettering	\$36.99
500-HMLAB-RED Red, "Infectious Linen" lettering	\$36.99
500-HMLAB-BLK Black, "Biohazardous Waste" lettering	\$36.99
500-HMLAB-GRN Green, "Resident Personal Linen" lettering	\$36.99

**HAMPER BAGS**

Item #	Price
500-HMBAG Nylon	\$38.99
500-HMBAG-A Antimicrobial	\$46.99
500-HMBAG-POLY-YEL Yellow poly disposable infectious linen	\$92.99
500-HMBAG-POLY-CLR Clear poly disposable	\$84.99
500-HMBAG-POLY-WHT White poly disposable soiled linen	\$95.99

**THREE PANEL PRIVACY SCREEN**

- 67" overall height
- Each panel measures 27" x 55"
- Three flame-retardant white vinyl panels
- White powder-coated steel frame
- WARRANTY: Limited

Item #	Price
<b>500-3PPS-A</b>	<b>\$96.95</b>
500-3PPS-C Privacy screen with 2" casters	\$114.95

**CHECK OUT**  
**promedsupply.com**

for full descriptions, color options and accessories

# Clinical & Diagnostic



## WELCH ALLYN® SPOT VITAL SIGNS®

- Measurements are quick, easy and accurate
- Lightweight and portable
- Large LCD display with backlight
- Includes 2 cuffs (1 adult size and 1 large adult size)
- Powered by a rechargeable 7.2 volt battery (included)
- **WARRANTY:** *Battery:* One year; *Spot vital signs unit:* Two years

Item #	Price
<b>WA42NTB-E1</b>	<b>\$2,039.99</b>
WA4700-60 Mobile stand	\$285.99
<b>WADS44-11</b>	<b>\$41.99</b>
Aneroid sphygmomanometer with adult one-piece cuff	
WADS44-12	\$43.99
Aneroid sphygmomanometer with adult one-piece cuff (large)	

**WA42NTB-E1**



**WelchAllyn®**

**WADS44-11**



**INNOVATIVE PRODUCT**

**500-ST61**



**500-ST154**



**TEMPERATURE  
MODULE**



**BLOOD  
PRESURE  
MODULE**



**PULSE OXIMETRY  
MODULE**



**CUSTOMER FAVORITE**

**SCHEDULE  
A DEMO**

Contact your territory manager or call us at **(800) 648-5190**

**FREE SHIPPING**  
ITEMS SHIP WITHIN 48 HOURS

**ADVIEW™ MODULAR  
DIAGNOSTIC STATION**

- Features individual blood pressure, temperature and pulse oximetry modules
- "Connect-and-go" capability between the blood pressure module and the temperature and pulse oximetry modules
- Automated and manual BP options
- Filac™ FasTemp™ provides 10-second oral or axillary temperatures
- SpO<sub>2</sub> module includes industry-leading Nellcor compatible technology
- **WARRANTY:** Two years

Item #	Price
<b>500-ST61</b> Diagnostic station	<b>\$1,699.99</b>
<b>500-ST154</b> Mobile stand	<b>\$285.99</b>

Item #	Price
500-ST62 Blood pressure module	<b>\$542.77</b>
500-ST63 Temperature and blood pressure modules	<b>\$956.96</b>
500-ST64 Pulse oximetry and blood pressure modules	<b>\$1,264.95</b>

600-SCHE



600-SWPP-2



600-SSPP



600-SPM

**Health o meter**<sup>®</sup>  
Professional

**CHECK OUT**  
[promedsupply.com](http://promedsupply.com)

for full descriptions on  
these and other items

**CUSTOMER FAVORITE**

**DIGITAL WHEELCHAIR  
RAMP SCALE**

- 1000 lb (454 kg) weight capacity
- Smart Load Cell technology
- Durable steel structure with non-slip platform
- WARRANTY: Three-year limited

Item #	Price
<b>600-SWPP-2</b>	<b>\$1,606.49</b>

**CHAIR  
SCALE**

- 440 lb (200 kg) weight capacity
- Graduation: 0.2 lb (0.1 kg)
- Seat: 18 1/4" x 14 1/2" x 17 1/2"
- Two wheels for ease of movement
- Flip arm and foot rest
- WARRANTY: Two-year limited

*Digital Only*

- EMR connectivity via USB
- 1 1/4" LCD display

Item #	Price
<b>600-SCHE</b>	<b>\$499.00</b>
Digital	
600-SCHM	<b>\$429.95</b>
Mechanical	

**BARIATRIC MEDICAL  
SCALE**

- 700 lb (310 kg) weight capacity
- Resolution: 0.2 lb (0.1 kg)
- Platform: 22" W x 15 3/4" D
- EMR connectivity via USB and serial ports
- 3/4" LCD display with 180° swivel for ease of viewing
- WARRANTY: Three-year limited

Item #	Price
<b>600-SSPP</b>	<b>\$1,007.99</b>

**PHYSICIAN'S  
BEAM SCALE**

- 500 lb (227 kg) weight capacity
- Large, easy-to-read measuring bar
- Non-slip, vinyl-covered platform
- WARRANTY: Two-year limited

Item #	Price
<b>600-SPM</b>	<b>\$183.00</b>



## INNOVATIVE PRODUCT

# A TRULY PORTABLE REAL-TIME ULTRASOUND BLADDER SCANNER



600-UBSU

Simple point-and-freeze operation makes the **Portascan+ LTC** the bladder scanner of choice for the busy clinical environment. The Portascan+ LTC calculates bladder volume automatically from an ultrasound image. A built-in thermal printer provides a permanent hard copy printout of results for accurate recordkeeping and archiving.

### PORTASCAN+ LTC BLADDER SCANNER

- Real-time ultrasound scanning and image capture
- Simple, intuitive software with touch-screen input
- 3.5-5.0 MHz abdominal probe for children and adults
- Automatic bladder volume calculation
- Touch-screen length and width calipers for kidney and other organ measurement
- Large database stores up to 250 images, records and reports
- Built in thermal printer for record keeping and archiving
- Requires a rechargeable battery (included) or external power supply
- WARRANTY: One year

Item #	Price
<b>600-UBSU</b>	<b>\$7,295.00</b>
Bladder scanner system and stand	
600-UBSU-P	\$2,458.00
Replacement convex 3.5-5.0 abdominal probe	
600-UBSU-BA	\$489.00
Rechargeable battery	
600-UBSU-GEL	\$12.99
Ultrasonic Gel	
600-UBSU-THP	\$79.99
Thermal paper	



## INTRODUCING A BREAKTHROUGH IN PT/INR TESTING

Coag-Sense PT/INR is the only point of care device that uses **direct detection of clot formation** to test rather than an electrochemical process.

The **Coag-Sense PT/INR** System provides accurate prothrombin time (PT) results in less than one minute, allowing abnormal results to be addressed immediately – no more losing valuable time waiting for lab results! Its **portable tabletop design** provides stability during testing and tests can be conducted using fresh capillary whole blood, venous or re-calcified plasma samples.



**600-PTINR-KIT**



### FREE METER WITH 5-BOX STARTER KIT

#### COAG-SENSE PT/INR MONITORING SYSTEM

- CLIA waived
- Direct clot detection technology for superb accuracy and precision
- Requires just 10 µ of blood
- Accurate and precise (CV 2.5%) results in one minute
- Measures the actual time required for clotting
- Barcoded test strips eliminate manual calibration
- Requires four AA batteries (included) or an A/C adapter (included)
- **WARRANTY:** One year with immediate replacement

Item #	Price
<b>600-PTINR-KIT</b>	<b>\$1,125.00</b>
New Customer Starter Kit includes: Meter Kit (600-PTINR-MTR) and 5 boxes of test strips (600-PTINR50)	
600-PTINR-MTR	\$1,445.00
Meter Kit includes: 1 meter, 1 user manual, 1 quick reference guide, 1 bag of auto single-usage lancets (quantity: 12), 4 AA alkaline batteries, 1 A/C power adapter, 1 minipipette two-pack, 1 rack of pipette tips (quantity: 96)	
600-PTINR-HL	\$65.00
5 high and low control strips, 2 control activation solution vials, pack of 14 sample transfer tubes	
600-PTINR-TT	\$10.50
PT/INR blood collection tubes (54/pk)	
600-PTINR50	\$289.00
PT/INR monitoring test strips (50/bx)	
040-970121	\$22.63
Haemolance Plus safety lancets (100/bx 48bx/cs)	

### FREE METER WITH 5-BOX STARTER KIT

- 1 meter (value: \$1,445.00)
- 5 boxes of test strips (quantity: 250; value: \$1,445.00)
- 10 high and low control strips
- 12 auto single-use lancets
- 1 minipipette two-pack
- 96 pipette tips

PRICE FOR  
STARTER KIT

**\$1,125<sup>00</sup>**

Total value of almost \$2,900!

# NEW CUSTOMER STARTER KITS



**WA04000-200**



**WA690**



**600-TH100**

## BRAUN THERMOSCAN® PRO 4000 EAR THERMOMETER

- ExacTemp™ technology ensures positioning and improved accuracy
- Provides fast results
- Automatic probe cover eject button
- Easy-to-read LCD display
- Requires two AA batteries (included)
- WARRANTY: Three-year limited

Item #	Price
<b>WA04000-200</b>	<b>\$189.75</b>
WA05074-800	\$17.27
Disposable probe covers (200/bx 4bx/cs)	

## NON-CONTACT INFRARED THERMOMETER

- Measure temperature by pointing thermometer at center of forehead
- Temperature measurement is quick, precise, hygienic and safe
- No probe covers necessary
- Requires four AAA batteries (included)
- WARRANTY: Two years

Item #	Price
<b>600-TH100</b>	<b>\$169.00</b>

## WELCH ALLYN® SURETEMP® PLUS 690 ELECTRONIC THERMOMETER

- Provides 4-second oral, 10-second rectal and 15-second adult axillary temperatures
- Monitor mode capability for measuring 3-5 minute temperatures
- Used probe covers can be ejected automatically
- Large LCD display
- Requires three AA batteries (not included)
- WARRANTY: *Instrument*: Two-year limited; *Probe*: One-year limited; *Probe well*: 90-day limited

Item #	Price
<b>WA690</b>	<b>\$272.99</b>
WA05031-105	\$11.39
Disposable probe covers (250/bx 20bx/cs)	

## BATTERIES

Item #	Price
500-BA-AA	\$0.42 (EA)
AA alkaline batteries	
500-BA-AAA	\$0.42 (EA)
AAA alkaline batteries	



## 3-BOX & 5-BOX SYSTEMS

The **CoaguChek XS PT/INR** system was developed to provide a high level of quality, accuracy and resident safety. Exclusive technology offers the flexibility and control that healthcare professionals need to confidently manage a wide range of residents on anticoagulation medication.

### Additional features:

- CLIA waived
- Accurate, precise results in about one minute
- Onboard controls – no external QC necessary
- Requires a small sample size (8 µ)
- Simple user interface with icon-based LCD
- Records 300 test results with date and time
- Auto-on with strip insertion; auto-off in three minutes
- Requires four AAA batteries (included)
- WARRANTY: One-year limited

To learn more, we invite you to contact your ProMed territory manager or give us a call at **(800) 648-5190**.

Item #	Description	Price
025-27001		<b>\$1,125.00*</b>
Includes: Analyzer, 3 boxes of test strips (025-15160), 1 box of CoaguChek Lancets		
025-56001		<b>\$1,325.00*</b>
Includes: Analyzer, 5 boxes of test strips (025-15160), 1 box of CoaguChek Lancets		
025-15160		<b>\$245.00</b>
Test strips with code chip (48/bx)		
025-50001		<b>\$73.53</b>
CoaguChek Lancets (200/bx)		

\*All promotional orders must be shipped directly to the new customer facility from Roche Diagnostics by June 30, 2013. Customers who have ordered the CoaguChek XS system product within the past 12 months will not be eligible for this promotion.



005-AS098-101



005-AS099-10



500-ASP1



**drive**  
MEDICAL RESPIRATORY & VENTILATION

**AIRSEP**



005-AS095-1

**CHECK OUT**  
**promedsupply.com**

for full descriptions on  
these and other items

## CUSTOMER FAVORITE

### VISIONAIRE OXYGEN CONCENTRATOR

- Compact and lightweight with near-silent operation
- Weight: 30 lbs
- Oxygen concentration: 1-5 LPM at 90% +5.5/-3%
- 20.8" H x 14.1" W x 11.5" D
- Operational temperature range: 41° F to 105° F
- WARRANTY: Parts: Five years

Item #	Price
<b>005-AS098-101</b>	<b>\$509.95</b>
005-AS098-104 w/Eco	\$549.95

### NEWLIFE INTENSITY 10 CONCENTRATOR

- Provides maximum pressure and flow
- Weight: 58 lbs
- Oxygen concentration: 2-9 LPM at 92% ±3%; 10LPM at 90% ±3%
- 27.5" H x 16.5" W x 14.5" D
- Operational temperature range: 41° F to 105° F
- WARRANTY: Five years

Item #	Price
<b>005-AS099-10</b>	<b>\$1,250.00</b>

### FREESTYLE BY AIRSEP

- The only truly wearable portable oxygen concentrator
- Weight: 4.4 lbs
- Oxygen concentration: 1-3 pulse settings; equivalent to continuous flow of 90% oxygen ±3%
- 8.6" H x 6.1" W x 3.6" D
- Operational temperature range: 41° F to 104° F
- WARRANTY: Unit: Three years; AirBelt: One year

Item #	Price
<b>005-AS095-1</b>	<b>\$3,395.00</b>

### HEAVY DUTY PORTABLE SUCTION MACHINE

- Controlled vacuum from 23" HG
- E-Z grip handle
- Base with cord wrap
- UL approved
- WARRANTY: Two-year limited

Item #	Price
<b>500-ASP1</b>	<b>\$148.95</b>

## CUSTOMER FAVORITE

### FINGERTIP PULSE OXIMETER (CLIP STYLE)



- Two-way display
- Not for continuous monitoring
- Rubber boot included
- Requires two AAA batteries (included)
- WARRANTY: Two-year limited

Item #	Price
500-PULSEOXF	\$49.95
500-BA-AAA AAA alkaline batteries	\$0.42 (EA)

### HANDHELD "HANDYOX" PULSE OXIMETER



- Customizable audible and visual alarm settings
- Individualized pulse rate and oxygen saturation high/low thresholds
- Durable rubber boot provides protection against accidental drops
- Requires two AA batteries (included)
- WARRANTY: Two-year limited

Item #	Price
500-PULSEOXH-DRV Includes finger probe	\$249.95
500-BA-AA AA alkaline batteries	\$0.42 (EA)

### PACIFICA ELITE NEBULIZER



- Weight: 4.1 lbs
- Ultra compact, portable design
- Powerful piston pump is efficient, dependable and quiet
- Complete with disposable neb kit
- WARRANTY: Five-year limited

Item #	Price
500-NC01M	\$19.99

### COMPARTMENT STYLE COMPRESSOR/NEBULIZER



- Maximum pressure: 30 psi
- Operating pressure: 12-14 psi
- Liter flow: 8 L/min
- 45° angle delivery allows resident to recline during treatment
- WARRANTY: Five-year limited

Item #	Price
500-NC01	\$39.95



# DISINFECT TO PROTECT

Surface disinfection plays a critical role in halting the spread of healthcare-associated infections (HAIs). Professional Medical stocks the items you need to protect both your residents and staff.

Whether you prefer sprays or wipes, our disinfectants meet EPA standards and are effective against a broad range of illness-causing microorganisms. We even stock products designed with special-needs populations in mind, including residents with respiratory sensitivities.

You can count on our infection prevention products to provide:


- The fast kill times you expect
- Compatibility with a wide variety of surfaces and equipment
- Ready-to-use formulations

To learn more, we invite you to contact your ProMed territory manager or give us a call at (800) 648-5190.

Item #	Description
600-SDWA-LG	Super Sani-Cloth Germicidal Disposable Wipe 6" x 6.75" (160/tub 12/cs)
600-SDWA-XL	Super Sani-Cloth Germicidal Disposable Wipe 7.5" x 15" (65/tub 6/cs)
600-SDWA-S	Super Sani-Cloth Germicidal Disposable Wipe Singles 11.5" x 11.75" (50/bx 3 bx/cs)
600-SDWAF-LG	Sani-Cloth AF3 Germicidal Disposable Wipe 6" x 6.75" (160/tub 12/cs)
600-SDWAF-XL	Sani-Cloth AF3 Germicidal Disposable Wipe 7.5" x 15" (65/tub 6/cs)
600-SDWAF-S	Sani-Cloth AF3 Germicidal Disposable Wipe Singles 11.5" x 11.75" (50/bx 3 bx/cs)
500-ESD03	Pure Hard Surface Cleaner 3oz (60/cs)
500-ESK11-Q6	Clearly Better Multi-Surface Cleaner 32oz spray (6/cs)

# Fall Prevention & Safety

## LANDING STRIP LOW PROFILE, BEVELED FLOOR MAT



- Injury-mitigating floor mat
- Stays in place even during care delivery
- Beveled edges accommodate wheelchairs and lifts
- Reduces tripping hazards

Item #	Price
600-FPB23-GR Gray, 23" x 68" x 1"	\$114.95
600-FPB23-BR Burgundy, 23" x 68" x 1"	\$114.95

## LANDING STRIP FLOOR MAT WITH NITE-GLO BORDER



- **Nite-Glo border** stands out in even dimly lit rooms to alert staff to the mat's presence
- Unusually thin and surprisingly firm
- Prevents tripping and permits carts, chairs, beds and tables to roll over easily
- Non-slip vinyl base

Item #	Price
600-FPB24BR-GLO 24" x 70" x 1"	\$139.95

## BEDSIDE SAFETY MAT



- One-piece design
- Helps reduce possible impact-related injuries
- Provides full coverage of bed exit area
- Made from high-impact foam
- Fire-resistant, wipe-clean cover

Item #	Price
600-FP55 24" x 66" x 2"	\$45.99

## BI-FOLD BEDSIDE SAFETY MAT



- Bi-fold design for easy storage
- Helps reduce possible impact-related injuries
- Provides full coverage of bed exit area
- Made from high-impact foam
- Skid-resistant bottom
- Fire-resistant, wipe-clean cover

Item #	Price
600-FP55-BF 24" x 66" x 2"	\$59.99

## CUSTOMER FAVORITE

### STRING ALARM



- Durable materials withstand institutional use
- Alarm settings inside battery compartment
- Unique alarm tone
- Unbreakable stainless steel mounting clip
- Requires one 9V battery (not included)
- **WARRANTY:** One year

Item #	Price
600-AA1	\$15.95
500-BA-9V 9 volt alkaline battery	\$1.17 (EA)

### BASIC ALARM



- Automatically resets when pressure is applied to sensor pad
- Off/low/high switch
- 0/2 second delay located in battery compartment
- Strain relief for sensor pad cable
- Compatible with our universal chair and bed sensors and seat belt systems
- Requires one 9V battery (not included)
- **WARRANTY:** One year

Item #	Price
600-AL1	\$29.95
500-BA-9V 9 volt alkaline battery	\$1.17 (EA)

### FLEX/ADVANCED ALARM



- Attaches to nurse call system
- Alarm is always in "on" position
- Automatically resets when pressure is applied to sensor pad
- Alarm settings inside battery compartment
- Compatible with our universal chair and bed sensors, seat belt systems and floor sensor pad (600-ALFP-36.5)
- Requires one 9V battery (not included)
- **WARRANTY:** One year

Item #	Price
600-AL2	\$49.95
500-BA-9V 9 volt alkaline battery	\$1.17 (EA)

### 3-IN-1 ULTRA ALARM



- Alarm is always in "on" position
- Automatically resets when pressure is applied to sensor pad
- Alarm settings inside battery compartment
- Nurse call compatible
- Compatible with our universal chair and bed sensors, seat belt systems and floor sensor pad (600-ALFP-36.5)
- Requires one 9V battery (not included)
- **WARRANTY:** One year

Item #	Price
600-AL3	\$54.95
500-BA-9V 9 volt alkaline battery	\$1.17 (EA)

**CUSTOMER FAVORITE**

**UNIVERSAL 45-DAY SENSOR PAD**



Sensor pads available for both chairs and beds.

- Works with many brands of fall monitors
- WARRANTY: 45 days

**Bed only:**

- Oversized to reduce false alarms
- 68" bed sensor cable

Item #	Price
600-ALCP-45 Chair sensor pad; 10" W x 15" L	\$17.95
600-ALBP-45 Bed sensor pad; 10" W x 28" L	\$20.20

**UNIVERSAL ONE-YEAR SENSOR PAD**



Sensor pads available for both chairs and beds.

- Moisture resistant
- Works with many brands of fall monitors
- WARRANTY: One year

**Bed only:**

- Oversized to reduce false alarms
- 68" bed sensor cable

Item #	Price
600-ALCP-365 Chair sensor Pad; 10" W x 15" L	\$31.15
600-ALPB-365 Bed sensor Pad; 10" W x 28" L	\$35.76

**CUSTOMER FAVORITE**

**ALARM/SENSOR FALL PREVENTION SYSTEM**



The fall prevention system includes the **Basic Alarm (600-AL1)** and either a

- **Universal One-Year Chair Sensor Pad (600-ALCP-365)**
- or
- **Universal One-Year Bed Sensor Pad (600-ALPB-365)**

Item #	Price
600-AL365-CSYS Chair system	\$49.95
600-AL365-BSYS Bed system	\$59.95

**FLOOR SENSOR PAD**



- 22" W x 39" L
- Use at the bedside or doorway to alert caregivers of resident activity
- Heavy-duty and ideal for institutional use
- Use with 3-in-1 Ultra (600-AL3) or Flex/Advanced (600-AL2) alarms

Item #	Price
600-ALFP-365	\$71.95

**INFRARED BEDSIDE MONITOR**



- Infrared sensor detects movement and triggers alarm
- Doesn't require wires, cords or pads
- Nurse call capability
- Multiple installation methods
- Requires one 9V battery (not included) or A/C adapter
- WARRANTY: One year limited

Item #	Price
500-ALIRM Infrared bedside monitor	\$54.55
500-ALIRM-BB Bed bracket	\$18.96
500-BA-9V 9 volt alkaline battery	\$1.17 (EA)

**SEATBELT MONITORING SYSTEMS**



- Works with the Ultra, Flex/Advanced and Basic alarms
- Alarms as soon as resident attempts to leave chair
- Mounts to most wheelchairs using existing hardware
- Large red loop is easy to release (600-ALBV only)

Item #	Price
600-ALBV Hook-and-loop seatbelt monitoring system	\$22.55
600-ALBB Buckle seatbelt monitoring system	\$24.55

**SINGLE CALL CORD**



Single call cord available in 8' and 10' lengths

- Durable materials withstand institutional use
- Ergonomically designed pendant is easy to grasp
- 1/4" phone plug works with most call systems
- Latex free

Item #	Price
600-ANCS-8 8 feet in length	\$8.99
600-ANCS-10 10 feet in length	\$14.88

**GERIATRIC PAD CALL CORD**



- 8' L
- Large pad is easily activated with an entire hand, fist or elbow
- Durable materials withstand institutional use
- 1/4" phone plug works with most call systems
- Latex free

Item #	Price
600-ANCS-8G	\$24.99



© Professional Medical, Inc.  
We reserve the right to correct any errors  
that may occur within this brochure.

Compliance, Smart Choice,  
Professional Medical, Inc., ProMed, and  
We Make Improving Care Easier! are  
trademarks of Professional Medical, Inc.

Orders will be shipped in the order they are received.  
Special priced products available while supplies last.  
Promotional Pricing is valid through June 30, 2013,  
unless otherwise noted. Orders placed after the  
promotional expiration date will be sold at our everyday  
low price. Applicable sales tax will be applied.  
Our shipping/handling policy procedures apply.

700-Smart Choice Equipment Solutions\_Q2\_13 [02/13]



**Professional Medical, Inc.**  
A Tradition of Quality, Value, & Trusted Service Since 1968

1917 Garnet Ct. New Lenox, IL 60451 ■ 800.648.5190 ■ fax 866.726.7416 ■ [promedsupply.com](http://promedsupply.com)

**Compliance®**  
Improve Care. Save Time. Reduce Costs.



**We make improving care easier!®**