# Equipment and Furniture Solutions: Summer 2013



PROFESSIONAL MEDICAL, INC.

# Equipment and Furniture Solutions: Summer 2013

#### **Mobility & Lifting**

#### 3-9

- Commodes
- Geri Chairs
- Fall Safety Gliders
- Fall Safety Mobility Chairs

• Cushions & Lap Positioners

## • Lifts & Slings

- Walkers & Rollators
- Wheelchairs & Accessories

#### **Beds & Furnishings**

- Armchairs
- Bariatric Beds Bedside Cabinets

• Chests & Wardrobes

- High-Low Beds
- Trapezes

#### Mattresses

- Foam Mattresses
- Alternating Pressure Mattress Systems
- Low Air Loss Mattress Systems

#### Fall Prevention & Safety 20 - 21 Alarms • Fall Safety Mats

• Call Cords • Pad, Floor & Belt Sensors

#### **Bathroom Safety, Carts & Hampers** 22 - 25

- Privacy Screens
- PVC Linen & Distribution Carts
- PVC Hampers • PVC Hydration Carts
- PVC Shower Chairs
- Wire & Poly Hampers

Scales

• Thermometers

#### **Clinical & Diagnostic**

- Bladder Scanners
- EKG/ECG Devices
- PT/INR Meters

#### Respiratory

- Concentrators
- Nebulizers
- Portable Suction Machines • Pulse Oximeters

• Vital Signs & Diagnostic Stations

Equipment and Furniture Solutions: Summer 2013

10 - 17 • Headboards & Footboards • Safety Rails & Pads 18 - 19

26 - 33

34 - 35

# Mobility & Lifting



### **STANDARD WHEELCHAIRS**

- 250 lb (16"); 300 lb (18"); 350 lb (20") weight capacities
- Powder-coated silver vein steel frame
- Embossed vinyl upholstery
- Durable urethane tires on composite wheels
- WARRANTY: Three-year limited

#### FIXED ARM

item #	Price
500-WC16F 16" width & swing-away footres	\$122.95 .ts
500-WC16E 16″ width & elevating legrests	\$132.95
500-WC18F 18" width & swing-away footres	\$114.95 ts
500-WC18E 18″ width & elevating legrests	\$122.95

#### PADDED DESK ARM - DETACHABLE

Item #	Price
500-WC16DF	\$125.95
16" width & swing-away footrests	5
500-WC16DE	\$138.95
16" width & elevating legrests	
500-WC18DF	\$125.95
18" width & swing-away footrest	5
500-WC18DE	\$138.95
18" width & elevating legrests	
500-WC20DF	\$144.95
20" width & swing-away footrests	5
500-WC20DE	\$149.95
20" width & elevating legrests	<b>•</b> ••••••
500-WC22DF	\$264.10
22" width & swing-away footrests	+
500-WC22DE	\$269.25
22" width & elevating legrests	<i><b>4</b>207.20</i>
500-WC24DF	\$319.95
24" width & swing-away footrests	
500-WC24DE	\$329.95
24" width & elevating legrests	ψυ27.75

### **OVERSIZED WHEELCHAIRS**

- 450 lb weight capacity
- Carbon steel frame with triple-coated chrome
- Double-embossed vinyl upholstery
- Durable urethane tires on composite wheels
- WARRANTY: Three-year limited


#### FULL ARM - DETACHABLE

ltem #	Price
500-WC16DF-FA	\$125.95
16″ width & swing-away footres	sts
500-WC16DE-FA 16″ width & elevating legrests	\$138.95
500-WC18DF-FA	\$132.95
18″ width & swing-away footres	sts
500-WC18DE-FA 18″ width & elevating legrests	\$138.95
500-WC20DF-FA	\$154.95
20″ width & swing-away footres	sts

Item #	Price
500-WC20DE-FA 20" width & elevating legrests	\$159.95
500-WC22DF-FA 22" width & swing-away footres	\$264.10 its
500-WC22DE-FA 22" width & elevating legrests	\$269.25
500-WC24DF-FA 24" width & swing-away footres	\$319.95 its
500-WC24DE-FA 24" width & elevating legrests	\$329.95



# Mobility & Lifting





#### STAND-ASSIST LIFT

- 400 lb weight capacity
- Boom height: 26.2" to 81.1"
- Base width: 24.5" closed; 40.2" open
- 4.5" base clearance
- Easy-roll 4" casters
- Four sling attachment points for a variety of stand-up sling designs
- WARRANTY: Frame: Lifetime; Electronic components and actuator: One year

Item #	Price
500-LP350-S	\$2,495.99
400 lb capacity	
500-S35S	\$99.99
Standing sling; 440 lb cap	acity

#### BARIATRIC BATTERY-POWERED LIFT

- 600 lb weight capacity
- Boom height: 32" to 74.5"
- Base width: 25.75" closed; 45.75" open
- 2" base clearance
- Six-point cradle
- Extra-wide coated cradle with 360° rotation
- WARRANTY: Frame: Lifetime; Electronic components and actuator: One year

Item #	Price
500-LP600	\$2,039.95
600 lb capacity	

# BATTERY-POWERED

- 450 lb weight capacity
- Boom height: 27.5" to 76"
- Base width: 15" closed; 41" open
- Low 2" base clearance
- For use with two-point and fourpoint slings
- Fast, easy charge from AC outlet
- WARRANTY: Lifetime limited

ltem #	Price
500-LP450 450 lb capacity	\$1,385.99
500-LP450-L	\$1,579.48
With removable, recharge 450 lb capacity	able battery;

# FOUR POINT AND STANDING

Item #	Price
<b>500-S45L-S</b> 600 lb capacity; Large	\$105.29
500-S45M-S 600 lb capacity; Medium	\$105.29
500-S45L-MC Mesh with commode; 600 lb capacity; Large	\$89.50
500-S45M-MC Mesh with commode; 600 lb capacity; Medium	\$89.50

# Mobility & Lifting



We offer a comprehensive line of cushions and accessories designed to help protect against skin breakdown, improve resident comfort, correct positioning and reduce sliding. The innovative APEX designs feature dual-density foam layers for pressure

redistribution and a dual gel-pod for increased comfort and stability. These highquality items are proudly made in the USA.





- Wedge shape for improved posture and slide control
- Breathable four-way stretch cover promotes active healing
- Slip-resistant bottom with quickrelease straps
- Meets CA 117 standards

Item #	Price
600-WC1816WGF-EA 18" x 16" x 4" to 2"	\$69.95
600-WC2016WGF-EA 20" x 16" x 4" to 2"	\$89.95

#### APEX LEG ABDUCTOR (POMMEL) GEL-FOAM CUSHION



- Contoured foam pommel for comfortable abduction
- Firm foundation eliminates hammocking
- Slip-resistant bottom with quickrelease straps
- Meets CA 117 standards

 Item #
 Price

 600-WC1816P-EA
 \$79.50

 18" x 16" x 3"

# BUY FIVE GET ONE FREE

Qualify by buying 5 or more of these cushions and lap positioners in any size or any combination *Promo code #CUSHB5G1* 

GEL-FOAM CUSHION

• High-density foam and a

proof, fire-resistant

600-WC1816D-EA

18" x 16" x 2" 600-WC1816-EA

18" × 16" × 3" 600-WC2016-EA

20" x 16" x 3" 600-WC2218-EA

22" x 18" x 3" 600-WC2418-EA

24" x 18" x 3"

Item #

dual-chamber gel-pod provide

• Low-shear, antimicrobial, fluid-

• Slip-resistant bottom with Velcro

• Meets CA 117 standards

straps to keep cushion in place

\$22.95

\$26.95

\$35.95

\$36.95

\$38.95

superior pressure relief and comfort

#### APEX GEL-FOAM CUSHION

- High-performance fabric promotes healing via moisture evaporation and cooling at pressure points
- Dual-chamber gel-pod reduces gel migration and eliminates bottoming out
- Breathable four-way stretch removable cover is antimicrobial and mildew-resistant with an anti-slip bottom
- Easy-to-attach quick-release straps
- Meets CA 117 standards

ltem #	Price
<b>600-WC1816DGF-EA</b> 18" x 16" x 2"	\$59.95
600-WC2016DGF-EA 20" x 16" x 2"	\$79.95

#### WHEELCHAIR CUSHION



- Combines a foam base and chambered gel insert for enhanced comfort and pressure relief
- Low-shear, antimicrobial, fire-resistant and fluid-proof removable cover
- Slip-resistant bottom with Velcro straps to keep cushion in place
- Meets CA 117 standards

Item #	Price
500-WC02V 16" x 18" x 2"	\$16.95
500-WC03V 16" x 18" x 3"	\$21.55
500-WC04V 16" x 18" x 4"	\$24.35



Itom #

#### LAP POSITIONERS

- Recommended for residents needing upper torso support or an elevated resting spot for hands
- Assist with upper-body alignment and increase safety by gently reminding residents to stay seated
- Available in self-releasing, full arm and desk arm styles

nem #	Price
600-WCLH	\$31.95
Self-releasing; 18-22" x 20" x 4"	
600-WCLH-F	\$31.95
Full arm; 27" x 18" x 3"	
600-WCLH-D	\$31.95
Desk arm; 27" x 18" x 3"	

### HEEL SUSPENSION



- Super-soft antimicrobial peachskin cloth protects thin, fragile skin
- Hook and loop closures allow boot to be fastened in multiple positions
- Extended, cushioned design distributes heel weight over a larger area, reducing pressure and increasing blood flow

ltem #	Price
600-PHE3-M Medium	\$21.95
600-PHE3-L Large/Extra Large	\$26.95

### HEEL FLOTATION BOOT





- Soft velour fabric
- Adaptable hook and loop closures
- Water-based gel pack for hot or cold therapy

ltem #	Price
600-PHG3	\$23.95
One size fits most (average adult)	

# **INNOVATIVE PRODUCT**



#### LIFTWALKER

- 300 lb weight capacity
- Ideal for individuals who need help sitting or standing
- Easy-to-deploy retractable standassist bars
- Simple-to-lock wheels
- WARRANTY: Lifetime limited

Item #	Price
600-AFLW	\$129.95
300 lb capacity	

# 2-BUTTON FOLDING WALKER AND ACCESSORIES



- 350 lb weight capacity
- Adjustable height: 32"-39"
- Easy push-button mechanisms
- U-shaped frame for greater clearance
- WARRANTY: Lifetime limited

Item #	Price
500-W101 300 lb capacity	\$27.95
<b>500-W101-5WF</b> 5" fixed wheels	\$30.25
500-W1A-5WF 5" fixed wheels with rear co	\$9.78 (PR) ap glides
500-W1A-GC Cap glides	\$5.30 (PR)
500-W1A-GT Tennis ball glides	\$11.98 (PR)

## 



- 300 lb weight capacity
- 6" black non-marring wheels
- Includes carry pouch under seat
- Deluxe loop locks
- Ergonomic handles
- WARRANTY: Lifetime limited

Item #	Price
500-RLT1-BL Blue	\$109.92
500-RLT1-RD Red	\$109.92

# FOLDING STEEL COMMODE



- 350 lb weight capacity
- Folds flat for convenient storage and transportation
- Durable plastic snap-on seat and lid
- Includes a 7.5-quart commode bucket, cover and splash shield
- WARRANTY: Lifetime limited

ltem #	Price
<b>500-C101</b> 350 lb capacity; seat with lid	\$48.19
500-C101-XW 650 lb capacity; 16.5" seat with	<b>\$99.99</b> lid
500-C101-P 12-qt replacement pail and cover	\$5.50

# Mobility & Lifting





for full descriptions on these and other items

### **INNOVATIVE PRODUCTS**

#### DYN-ERGO FALL SAFETY SCOOT CHAIR

- Provides safer foot propulsion for those at risk of falling in a conventional wheelchair
- Enables correct positioning with customized adjustability
- Addresses mobility, comfort, dignity and positioning needs
- Non-institutional appearance
- WARRANTY: Extended

ltem #	Price
Specify vinyl upholstery	
<b>500-CPS18.5</b> 250 lb capacity	\$999.00
500-CPS24 400 lb capacity	\$1,693.00
500-CPS-SWFR Swing-away footrests	\$200.00 (PR)

#### EVOLUTION FALL SAFETY MOBILITY CHAIR

- Permits safe self-propulsion using one's feet and/or arms
- Helps maintain mobility and dignity
- Highly maneuverable
- and comfortable
- Non-institutional appearance
- WARRANTY: Extended

# Item # Price Specify vinyl upholstery \$ 500-RC18.5 \$1,471.00 250 lb capacity \$

#### COMP-POSITION FALL SAFETY CHAIR

- Designed for those who require support and positioning during periods of extended sitting
- Tilt action helps reduce pressure for better skin management
- Can be reconfigured for assignment to other users
- Non-institutional appearance
- WARRANTY: Extended

Item #	Price
Specify vinyl upholster	y
500-RTS	\$2,158.00
250 lb capacity	



those receiving memory care, to benefit from the therapeutic effects of a rocking chair while minimizing the risk of falls. A patented auto-lock mechanism permits safe transfers.



### FALL SAFETY GLIDERS

# THERA-GLIDE DECO FALL SAFETY GLIDER

ltem #	Price
Specify wood stain & ving	yl upholstery
<b>500-CGL20-DCB</b> 275 lb capacity; Classic back	\$919.00
500-CGL20-DCPB	<b>\$954.00</b>
275 lb capacity; Comfort Plus	back
500-CGL20-DSB	\$919.00
275 lb capacity; Seamless Simpli	city back
500-CGL24-DCB 450 lb capacity; Classic back	\$1,251.00
500-CGL24-DCPB	\$1,286.00
450 lb capacity; Comfort Plus	back
500-CGL24-DSB	\$1,251.00
450 lb capacity; Seamless Simpli	icity back

# THERA-GLIDE MISSION FALL SAFETY GLIDER

Item #	Price
Specify wood stain & vinyl upholstery	
500-CGL20-MCB 275 lb capacity; Classic back	\$884.00
<b>500-CGL20-MCPB</b>	\$919.00
275 lb capacity; Comfort Plus	back
500-CGL20-MSB	\$884.00
275 lb capacity; Seamless Simpli	city back
500-CGL24-MCB 450 lb capacity; Classic back	\$1,216.00
500-CGL24-MCPB	\$1,251.00
450 lb capacity; Comfort Plus	back
500-CGL24-MSB	\$1,216.00
450 lb capacity; Seamless Simpli	city back

# THERA-GLIDE UPHOLSTERED FALL SAFETY GLIDER

ltem #	Price
Specify wood stain & viny	l upholstery
500-CGL20-USCB 275 lb capacity; Classic back	\$884.00
500-CGL20-USCPB 275 lb capacity; Comfort Plus	\$919.00 back
500-CGL20-USSB\$884.00275 lb capacity; Seamless Simplicity back	
500-CGL24-USCB 450 lb capacity; Classic back	\$1,216.00
500-CGL24-USCPB 450 lb capacity; Comfort Plus	\$1,251.00 back
500-CGL24-USSB 450 lb capacity; Seamless Simplie	\$1,216.00 city back

# THERA-GLIDE TRANSITIONAL FALL SAFETY GLIDER

ltem #	Price
Specify wood stain & vinyl upholstery	
500-CGLT23	\$1,199.00
350 lb capacity	

• Compatible with all styles of the Thera-Glide Fall Safety Gliders

Item #	Price
Specify vinyl upholstery	
500-CGLRA-28	\$99.00
Undercarriage cover	

# Beds & Furnishings



# When it's time to renovate your facility, turn to the experts.

Professional Medical has partnered with Basic American to create a complimentary design program that aims to enhance resident comfort and wellness. We want to help you create a warm, home-like environment that your residents and their families will enjoy.

Not only will you avoid having to pay an interior designer's hourly charges, you'll also be partnering with a team of professionals who are in tune with the specific challenges of designing for long-term care. Our team is detail-oriented and aware of the codes and regulations that need to be taken into account when designing a space that is both attractive and functional. We provide the following services, all at no charge:

- Space planning
- Schematic design
- Design development
- Furniture selection
- Color and materials selection
- Custom furniture and millwork design
- Artwork selection
- Graphics and signage
- Accessories and lighting specifications
- Installation administration
- In-house purchasing



Whether you're looking to spruce up or completely overhaul resident rooms, dining areas, activity spaces, common areas or office settings, our experienced interior designers are ready to help you create the perfect look for your facility.

To learn more, contact your ProMed territory manager or call our Customer Care team at (800) 648-5190.

# Beds & Furnishings



#### LIBERTY BED

- 450 lb weight capacity
- 76" or 80" L x 35" W
- Adjustable brackets allow bed to be used at both 76" and 80" lengths
- Height range: 9.9" to 24.4"
- 360-degree mobility at any height
- Quiet, smooth operation with three DC motors
- Surface treated with antimicrobial SilverSolutions
- WARRANTY: Three years

Item #	Price
075-LBT7174	\$735.95
Liberty Bed (assist rails sold sep	parately)
075-AB53600	<b>\$59.95</b>
Headboard & footboard; Sola	r Oak
500-BSC3-O	<b>\$415.99</b>
3-drawer bedside cabinet; Sol	ar Oak
500-DRS3-O 4-drawer chest; Solar Oak	\$523.99

#### BARIATRIC



### FREE 42" FOAM MATTRESS AND PAIR OF "T" RAILS

#### FULL ELECTRIC BARIATRIC BED

- 600 lb weight capacity
- 88" L x 42" W
- Height range: 18.25" to 26.25"
- Heavy-duty frame for strength and resident safety
- Split pan design with removable bed ends for easy setup

# **CUSTOMER FAVORITE**

- Quiet, smooth operation
- WARRANTY: Frame: Five years; Other parts and components: One year; Welds: Lifetime

Item #	Price
500-BB600-PKG	\$2,249.00
Package	



#### BARIATRIC EXPANDABLE BED

- 750 lb weight capacity
- 80"-88" L × 39"-48" W
- Height range: 15.75" to 29.25"
- Smokey Hollow headboard and footboard available in Solar Oak laminate
- WARRANTY: Frame: 15 years; Electronic components: Two years; Welds: Lifetime

Item #	Price
600-BED4 750 lb capacity	\$3,149.65
600-BED4-HB1	\$77.50
Headboard; Rich Solar Oak	laminate
600-BED4-FB1	\$77.50
Footboard; Rich Solar Oak I	aminate

#### ZENITH 9000 BED

- 600 lb weight capacity
- 76" or 80" L x 35" W
- Adjustable brackets allow bed to be used at both 76" and 80" lengths
- Height range: 7" to 30"
- 360-degree mobility at any height
- Fast, even-speed Timesaver Technology™ motors
- Surface treated with antimicrobial SilverSolutions
- WARRANTY: Frame: 15 years; Motors/control box: Five years; Other electronics: Three years

Item #	Price
075-ZZ98174 Assist rails sold separately	\$1,285.95
075-Z04500CCR Headboard & footboard; Wild	\$61.25 I Cherry
500-BSC3-C 3-drawer bedside cabinet; W	\$415.99 ild Cherry
500-DRS3-C 4-drawer chest; Wild Cherry	\$523.99



### **INNOVATIVE PRODUCT**



#### ZENITH SLIDE-W-I-D-E

- 600 lb weight capacity
- Quickly and easily adjusts to 35", 39" and 42" widths
- Ideal for bariatric residents and those accustomed to wider surfaces
- Fast adjustment to fit through narrow doorways
- WARRANTY: Frame: 15 years; Motors and control box: 5 years

ltem #	Price
075-ZZ98674 600 lb capacity	\$1,599.95
075-ZW48500 Headboard & footboard;	\$85.25 Medium Oak





#### **HIGH-LOW BEDS**

#### BASIC HIGH-LOW BED

- 450 lb weight capacity
- 80" L x 35" W
- Height range: 8" to 26"
- Roll-at-any-height design
- Locking system
- WARRANTY: Frame: 15 years; Electronic components: Two years; Welds: Lifetime

Item #	Price
600-BED1	\$998.65
450 lb capacity	

#### BASIC<sup>+</sup> HIGH-LOW BED

- 450 lb weight capacity
- 80" L x 35" W
- Height range: 8" to 26"
- Advanced caster locking system
- Optional 42" expandable deck
- WARRANTY: Frame: 15 years; Electronic components: Two years; Welds: Lifetime

ltem #	Price
600-BED2	\$1,138.65
450 lb capacity	

#### ULTIMATE HIGH-LOW BED

- 500 lb weight capacity
- Expands from 76" to 84" L and 35" to 42" W
- Height range: 6.88" to 30.88"
- Roll-at-any-height design
- Provides more workable surface for repositioning and care
- WARRANTY: Frame: 15 years; Electronic components: Two years; Welds: Lifetime

 Item #
 Price

 600-BED3
 \$1,618.65

 500 lb capacity
 \$1

#### SMOKEY HOLLOW HEADBOARDS AND FOOTBOARDS

ltem #	Price
600-BED-HB1 Headboard; Solar Oak laminate	\$42.60
600-BED-FB1 Footboard; Solar Oak laminate	\$42.60
600-BED-HB2 Headboard; Amber Cherry lamin	\$42.60 ate
600-BED-FB2 Footboard; Amber Cherry lamina	\$42.60 te

#### MILL CREEK HEADBOARDS AND FOOTBOARDS

Item #	Price
600-BED-HB3 Headboard; Solar Oak laminate	\$42.60
600-BED-FB3 Footboard; Solar Oak laminate	\$42.60
600-BED-HB4 Headboard; Amber Cherry lamir	\$42.60 nate
600-BED-FB4 Footboard; Amber Cherry lamina	\$42.60 te

# AUTO-TOUCH OVERBED TABLE • Adjustable from 28" to 45" • 15" x 30" tabletop • Walnut wood grain laminated top • 2" swivel casters for easy maneuverability • WARRANTY: Lifetime limited \$49.99 500-TTBED 15" x 30" tabletop HALF AND FULL VINYL BED RAIL PADS TRAPEZE BAR • 34" or 69" L x 11" W x 1" D • Half-length and full-length styles available • Provide protective padding to reduce the risk of resident injury Price 600-BRP1134-PR \$21.25 Half-length vinyl bed rail pads; rails sold separately \$37.95

600-BRP1169-PR Full-length vinyl bed rail pads; rails sold separately

# SAFETY HALF-LENGTH RAILS



- Fits all spring-style healthcare beds
- No-gap style
- Steel construction
- WARRANTY: Lifetime limited

	Item #	Price
)	500-RAIL-H Half-length rails	\$56.95

#### BARIATRIC



- 750 lb weight capacity
- Heavy-duty grab handle and multiheight position strap
- Easy-to-use rung-style strap attachments for quick relocation
- Freestanding design

Price 600-T01B-F \$395.55 Freestanding

# New Equipment & Furniture Now

Tailored leasing plans to fit your needs and budget

Professional Medical's exclusive Loan-To-Own (LTO) **Program** can bolster your facility's therapeutic equipment inventory, help you acquire new state-of-the-art merchandise and save you thousands of dollars over other traditional rental programs.

#### Use your funds wisely now, not during survey time!

Allow Professional Medical to help evaluate the equipment needed to provide the safety your employees expect and the quality care your residents deserve. Now, with our LTO solution, you can afford the equipment needed to provide top-notch care while remaining compliant with state and federal guidelines.

Your ProMed sales representative can walk you through preparing a customized budget analysis offering the best therapeutic equipment backed by the most cost-effective leasing program in the industry.



# Beds & Furnishings

500-BSC1-C



# BEDSIDE CABINETS

CHERRY

OKLAHOMA OAK

WHITE ROCK MAPLE

- 29" W x 18" D x 29.5" H
- Durable nickel-finish hardware
- All rounded corners for safety
- WARRANTY: Five-year limited

Item #	Price
<b>500-BSC1-C</b>	\$129.99
1-door, 1-drawer cabinet; Cherry	/
500-BSC1-O	<b>\$129.99</b>
1-door, 1-drawer cabinet; Oklah	oma Oak
500-BSC1-M	\$129.99
1-door, 1-drawer cabinet; White	Rock Maple
<b>500-BSC2-C</b> 3-drawer cabinet; Cherry	\$138.99
500-BSC2-O	\$138.99
3-drawer cabinet; Oklahoma Oc	1k
500-BSC2-M	\$138.99
3-drawer cabinet; White Rock M	aple

# 3-DRAWER AND 4-DRAWER CHESTS

- 29" W x 18" D x 29.5" H (3-dr.); 29" W x 18" D x 37" H (4-dr.)
- Durable nickel-finish hardware
- All rounded corners for safety
- WARRANTY: Five-year limited

ltem #	Price
500-DRS1-C 3-drawer chest; Cherry	\$211.99
500-DRS1-O 3-drawer chest; Oklahoma Oak	\$211.99
500-DRS1-M 3-drawer chest; Maple	\$211.99
<b>500-DRS2-C</b> 4-drawer chest; Cherry	\$222.99
500-DRS2-O 4-drawer chest; Oklahoma Oak	\$222.99
500-DRS2-M 4-drawer chest; White Rock Map	\$222.99 Ile

# 1-DOOR WARDROBES

- 22" W x 23" D x 78" H
- Durable nickel-finish hardware
- All rounded corners for safety
- WARRANTY: Five-year limited

Item #	Price
500-WRD1-2-C	\$386.99
1-door, 2-drawer wardrobe; C	Cherry
500-WRD1-2-O 1-door, 2-drawer wardrobe; C	\$386.99 Dklahoma Oak
500-WRD1-2-M 1-door, 2-drawer wardrobe; W	\$386.99 hite Rock Maple

# 2-DOOR WARDROBES

500-WRD2-C

- 30" W x 23" D x 78" H
- Durable nickel-finish hardware
- All rounded corners for safety
- WARRANTY: Five-year limited

ltem #	Price
<b>500-WRD2-C</b> 2-door wardrobe; Cherry	\$442.99
500-WRD2-O	\$442.99
2-door wardrobe; Oklahoma Oc	<sup>ak</sup>
500-WRD2-M	\$442.99
2-door wardrobe; White Rock N	Naple
500-WRD2-2-C	\$484.99
2-door, 2-drawer wardrobe; Che	erry
500-WRD2-2-O	<b>\$484.99</b>
2-door, 2-drawer wardrobe; Okl	ahoma Oak
500-WRD2-2-M	\$484.99
2-door, 2-drawer wardrobe; White	Rock Maple



FINISH AND	FABRICS AS S	HOWN						
FINISH	CRYPTON F	ABRIC - CURL	EY Q					
Natu Mah	ral ogany	27	er er	C	<u>-</u> 2)	2)	C	
	Incense	Henna	laurel	Auburn	lark	Sage	Granite	

All Kwalu styles are available in a wide variety of finishes and fabric options.

Please contact us for more details.

#### TIVOLI ARMCHAIR

- Curved back, custom carved arms and split-back detailing
- Superior cushioning
- Easy to clean
- WARRANTY: Construction and finish: 10 years

Item #	Price
500-AC251111	\$499.88
Tivoli	

#### SYMPHONY ARMCHAIR

- Versatile and appropriate for a wide range of environments
- Scuff-, scratch- and water-resistant finish with antimicrobial product protection
- Arms engineered for easy ingress and egress
  WARRANTY: Construction and
  - finish: 10 years

```
        Item #
        Price

        500-AC32169
        $454.88

        Symphony
        $454.88
```

#### REGAL ALBERT ARMCHAIR

- Traditional styling
- Scuff-, scratch- and water-resistant finish with antimicrobial product protection
- WARRANTY: Construction and finish: 10 years

ltem #	Price
500-AC54124	\$429.88
Regal Albert	

#### BEDSIDE CABINETS

- Softly rounded edges and waterfall detailing
- Inset doors and drawers
- Ideal for all healthcare interiors
- WARRANTY: Construction and finish: 10 years

Item #	Price
500-BSC1-S11R 1-door, 1-drawer; Mahogany	\$249.88
500-BSC3-S30	\$269.88
3-drawer; Mahogany	


#### 600-M700-76

# 3-YEAR PRESSURE REDISTRIBUTION MATTRESS

- 300 lb weight capacity
- Three pressure zones
- Channel-cut high-density polyurethane foam
- WARRANTY: Three years, non-prorated

Item #	Price
<b>600-M700-76</b> 36" x 76" x 6"	\$139.95
600-M700-80 36" x 80" x 6"	\$139.95
600-M700-84 36" x 84" x 6"	\$149.95

# DYNAMIC ELITE MATTRESS WITH RAISED RAILS



- 400 lb weight capacity
- Quickly channels heat away from the body, helping to reduce pressure ulcer incidence
- Provides superior pressure redistribution
- Side rails decrease fall risk • WARRANTY: 10 years, non-prorated

1	
Item #	Price
600-M1000-76RR 36" x 76" x 6"	\$255.99
600-M1000-80RR 36" x 80" x 6"	\$264.99
600-M1000-84RR 36" x 84" x 6"	\$269.99



# 6-YEAR PRESSURE REDISTRIBUTION MATTRESS

- 350 lb weight capacity
- 4-way stretch nylon waterproof cover
- Provides additional pressure reduction in heel area
- Solid-edge support for easy resident transfer
- WARRANTY: Six years

ltem #	Price
600-M850-76 36″ x 76″ x 6″	\$187.95
600-M850-80 36" x 80" x 6"	\$187.95
600-M850-84 36" x 84" x 6"	\$198.95

#### **CUSTOMER FAVORITE**

15-YEAR PRESSURE REDISTRIBUTION MATTRESS



- 450 lb weight capacity
- 4-way stretch nylon waterproof cover
- Designated head and foot sections
- Concealed zipper with Barrier Stop Över-Flap
- WARRANTY: 15 years

ltem #	Price
600-M900-76 36″ x 76″ x 6″	\$214.99
600-M900-80 36" x 80" x 6"	\$214.99
600-M900-84 36" x 84" x 6"	\$225.99

• 400 lb weight capacity • Quickly channels heat away

BUY

Promo code #FOAMB5G1

DYNAMIC ELITE MATTRESS

GET ONE FREE

more foam mattresses of

- from the body, helping to reduce pressure ulcer incidence
- Provides superior pressure redistribution
- WARRANTY: 10 years, non-prorated

Item #	Price
600-M1000-76 36" x 76" x 6"	\$229.98
600-M1000-80 36" x 80" x 6"	\$229.98
600-M1000-84 36" x 84" x 6"	\$245.99

#### **CUSTOMER FAVORITE**

# ALTERNATING PRESSURE MATTRESS SYSTEM

- 350 lb weight capacity
- Includes 20-bladder mattress and high-flow pump
- Helps prevent and treat pressure sores
- Adjustable comfort settings allow pressure to be customized by resident weight
- WARRANTY: Pump: One year

Item #	Price
600-AM7000 36" x 80" x 8"	\$465.95

#### TRUE LOW AIR LOSS MATTRESS SYSTEM WITH PULSATE



- 1000 lb weight capacity
- Includes mattress and LAL/ pulsate pump
- Provides a cool and dry sleeping surface
- Exclusive pulsate mode reduces pressure and increases blood flow
- WARRANTY: Limited

ltem #	Price
600-AM9000	\$2,015.99
36" x 80" x 10"	



- 350 lb weight capacity
- Provides low air loss and alternating pressure therapy in 10,15, 20 and 25-minute cycles
- 10 adjustable comfort cycles allow pressure to be customized by resident weight
- Pillow top keeps head and neck stable
- WARRANTY: Limited

Item #	Price
600-AM8000	\$789.95
36" x 80" x 8"	

#### BARIATRIC LOW AIR LOSS WITH ALTERNATING PRESSURE MATTRESS



- 1000 lb weight capacity
- Provides low air loss and alternating pressure therapies
- Microprocessor unit provides continuous feedback to equalize pressure and maximize comfort
- WARRANTY: Limited

Item #	Price
600-AM9000B 42" x 80" x 10"	\$2,189.99
600-AM9000B48 48" x 80" x 10"	\$2,209.99



# DISINFECT TO PROTECT

Surface disinfection plays a critical role in halting the spread of healthcare-associated infections (HAIs). Professional Medical stocks the items you need to protect both your residents and staff.

Whether you prefer sprays or wipes, our disinfectants meet EPA standards and are effective against a broad range of illness-causing microorganisms. We even stock products designed with special-needs populations in mind, including residents with respiratory sensitivities.


#### You can count on our infection prevention products to provide:

- The fast kill times you expect
- Compatibility with a wide variety of surfaces and equipment
- Ready-to-use formulations

# To learn more, we invite you to contact your ProMed territory manager or give us a call at (800) 648-5190.

Item #	Descrption
600-SDWA-LG	Super Sani-Cloth Germicidal Disposable Wipe 6" x 6.75" (160/tub 12/cs)
600-SDWA-XL	Super Sani-Cloth Germicidal Disposable Wipe 7.5" x 15" (65/tub 6/cs)
600-SDWA-S	Super Sani-Cloth Germicidal Disposable Wipe Singles 11.5" x 11.75" (50/bx 3 bx/cs)
600-SDWAF-LG	Sani-Cloth AF3 Germicidal Disposable Wipe 6" x 6.75" (160/tub 12/cs)
600-SDWAF-XL	Sani-Cloth AF3 Germicidal Disposable Wipe 7.5" x 15" (65/tub 6/cs)
600-SDWAF-S Sani-Cloth AF3 Germicidal Disposable Wipe Singles 11.5" x 11.75" (50/bx 3 bx/cs)	
500-ESD03	Pure Hard Surface Cleaner 3oz (60/cs)
500-ESK11-Q6	Clearly Better Multi-Surface Cleaner 32oz spray (6/cs)

# Fall Prevention & Safety



#### UNIVERSAL 45-DAY SENSOR PAD



Sensor pads available for both chairs and beds.

- Works with many brands of fall monitors
- WARRANTY: 45 days

#### Bed only:

- Oversized to reduce false alarms
- 68" bed sensor cable

ltem #	Price
600-ALCP-45	\$17.95
Chair sensor pad; 10" W x 15" L	
600-ALPBP-45	\$20.20
Bed sensor pad; 10" W x 28" L	



Sensor pads available for both chairs and beds.

- Moisture resistant
- Works with many brands of fall monitors
- WARRANTY: One year

#### Bed only:

- Oversized to reduce false alarms
- 68" bed sensor cable

ltem #	Price
600-ALCP-365	\$31.15
Chair sensor pad; 10" W x 15" L	
600-ALPB-365	\$35.76
Bed sensor pad; 10" W x 28" L	



The fall prevention system includes the **Basic Alarm** (600-AL1) and either a

- Universal One-Year Chair Sensor Pad (600-ALCP-365) or
- Universal One-Year Bed Sensor Pad (600-ALPB-365)



FLOOR SENSOR PAD

- 22" W x 39" L
- Use at the bedside or doorway to alert caregivers of resident activity
- Heavy-duty and ideal for institutional use
- Use with 3-in-1 Ultra (600-AL3) or Flex/Advanced (600-AL2) alarms

ltem #	Price
600-AL365-CSYS	\$49.95
Chair system	
600- AL365-BSYS	\$59.95
Bed system	

Item #	Price
600-ALFP-365	\$71.95
22" W x 39" L	

# INFRARED BEDSIDE



- Infrared sensor detects movement and triggers alarm
- Doesn't require wires, cords or pads
- Nurse call capability
- Multiple installation methods
- Requires one 9V battery (not included) or A/C adapter
- WARRANTY: One-year limited

Item #	Price
500-ALIRM Infrared bedside monitor	\$54.55
500-ALIRM-BB Bed bracket	\$18.96
500-BA-9V 9 volt alkaline battery	\$1.17 (EA)

SEATBELT MONITORING SYSTEMS
4



- Works with the Ultra, Flex/ Advanced and Basic alarms
- Alarms as soon as resident attempts to leave chair
- Mounts to most wheelchairs using existing hardware
- Large red loop is easy to release (600-ALBV only)

Item #	Price
600-ALBV	\$22.55
Hook and loop seatbelt monitor	ring system
600-ALBB	\$24.55
Buckle seatbelt monitoring syste	em

# SINGLE CALL



Single call cord available in 8' and 10' lengths

- Durable materials withstand institutional use
- Ergonomically designed pendant is easy to grasp
- ¼" phone plug works with most call systems
- Latex free

Item #	Price
600-ANCS-8 8 feet in length	\$8.99
600-ANCS-10 10 feet in length	\$14.88

#### GERIATRIC PAD CALL CORD



- 8' length
- Large pad is easily activated with an entire hand, fist or elbow
- Durable materials withstand institutional use
- ¼" phone plug works with most call systems
- Latex free

Item #	Price
600-ANCS-8G	\$24.99
8 feet in length	



#### **CUSTOMER FAVORITE**

# Bathroom Safety, Carts & Hampers



Our shower chairs and hydration/ distribution carts are made from durable, heathcare-grade PVC that inhibits the growth of bacteria and is free of phthalates and lead.

These products are reinforced at all stress-related areas for increased durability and performance. They won't rust, fade or stain!

#### PVC SHOWER CHAIR WITH CASTERS

- 300 lb weight capacity
  Internal width: 18";
- external width: 22"3" or 5" casters
- Push-pull handle for easy

5" casters

- maneuvering
- Fast-drying removable mesh slingWARRANTY: PVC frame: Five-year
- limited
  Item # Price
  Specify mesh color
  500-PVC-SC18 \$99.95
  3" casters
  500-PVC-SC18-C5 \$136.95

#### PVC SHOWER CHAIR WITH FOLDING FOOTREST

- 300 lb weight capacity
- Internal width: 18"; external width: 22"
- 3" twin rust-proof casters
- Folding footrest
- Fast-drying removable mesh sling
- WARRANTY: PVC frame: Five-year limited

Price
\$135.95

#### PVC BARIATRIC SHOWER CHAIR

- 425 lb weight capacity
- Internal width: 26"; external width: 30"
- 4" rust-proof casters
- Push-pull handle for easy maneuvering
- Fast-drying removable mesh sling
- WARRANTY: *PVC frame:* Five-year limited

Item #	Price
Specify mesh color	
500-PVC-SC26	\$195.95
425 lb capacity	



#### 36-QUART HYDRATION CART

- 36-quart ice chest
- Ice scoop
- Storage shelf
- 3" twin rust-proof casters
- WARRANTY: *PVC frame:* Five-year limited

Item #	Price
Specify mesh or solid vinyl color	
500-PVC-HC	\$192.99
36 qt capacity	

#### 48-QUART HYDRATION CART

- 48-quart ice chest
- Ice scoop
- Storage shelf
- Pull-out storage/drain tub
- 3" twin rust-proof casters
- WARRANTY: *PVC frame:* Five-year limited

Item #	Price
Specify mesh or solid vinyl color	
500-PVC-HC2	\$214.99
48 qt capacity	

Our **fixed height bath** is ideal for longterm care use, especially when transfer from a wheelchair is required. The powered traverse makes it easy for residents to get in and out of the bath, helping them retain their independence. The seat can be easily lowered outside the tub to a convenient height for comfort and safety.

The bath also features high-end touches such as a stylish, extra-deep tub and a convenient hand control. Battery back-up comes standard for total peace of mind.

#### **Additional features:**

- 330 lb lifting capacity
- 60" and 66" lengths available
- Left- and right-hand models available
- Seat can be set at any height
- Smooth, easy-clean surface
  - WARRANTY: Two years

 Item #
 Price

 500-TUB1
 \$10,497.00

 330 lb lifting capacity
 \$10,497.00



# THE FIXED HEIGHT BATH

# Bathroom Safety, Carts & Hampers



#### HAMPER WITH MESH ZIPPER BAG

- Zipper-front opening eliminates lifting heavy loads over the top of the hamper
- Durable high-impact plastic lid
- 3" deluxe, heavy-duty reinforced casters
- Floor plate to support retainer bag
- WARRANTY: PVC frame: Five-year limited

Item #	Price
Specify mesh or solid vinyl color	
500-PVC-HM55-1 Single hamper; 26 gallon	\$129.99
500-PVC-HM55-2 Double hamper; 52-gallon	\$210.99

#### SOLID SHELF LINEN CART

- 100 lb carrying capacity per shelf
- Shelves are 12" H x 20" W (see below for available lengths)
- 10" W handles
- WARRANTY: *PVC frame:* Five-year limited

ltem #	Price
Specify mesh or solid vin	yl color
500-PVC-LC25-2 2-shelf cart; 20" x 25"	\$154.99
500-PVC-LC32-2 2-shelf cart; 20" x 32"	\$194.86
500-PVC-LC45-2 2-shelf cart; 20" x 45"	\$252.99
500-PVC-LC50-2 2-shelf cart; 20" x 50"	\$268.99
<b>500-PVC-LC25-3</b> 3-shelf cart; 20" x 25"	\$183.72
500-PVC-LC32-3 3-shelf cart; 20" x 32"	\$254.99
500-PVC-LC45-3 3-shelf cart; 20" x 45"	\$278.99
500-PVC-LC50-3 3-shelf cart; 20" x 50"	\$342.99

#### PERFORATED SHELF LINEN CART

- 75 lb carrying capacity per shelf
  Shelves are 12" H x 20" W (see
- below for available lengths)
- 10" W handles
- WARRANTY: *PVC frame:* Five-year limited

ltem #	Price
Specify mesh or solid vin	yl color
500-PVC-LC25-2A 2-shelf cart; 20" x 25"	\$154.99
500-PVC-LC32-2A 2-shelf cart; 20" x 32"	\$194.86
500-PVC-LC45-2A 2-shelf cart; 20" x 45"	\$243.58
500-PVC-LC50-2A 2-shelf cart; 20" x 50"	\$277.84
<b>500-PVC-LC25-3A</b> 3-shelf cart; 20" x 25"	\$200.38
500-PVC-LC32-3A 3-shelf cart; 20" x 32"	\$242.04
500-PVC-LC45-3A 3-shelf cart; 20" x 45"	\$340.21
500-PVC-LC50-3A 3-shelf cart; 20" x 50"	\$342.27

# DISTRIBUTION

- Shelf: 20" W x 43.63" L
- Length with handles: 50"; height with casters: 66.38"
- Hanging rack: 36.75"
- Bin:16" W x 12.25" L x 5.5" D
- WARRANTY: *PVC frame*: Five-year limited

Specify mesh or solid vinyl color	
500-PVC-LC67-6A	\$334.99
Distribution cart	



# SINGLE "EASY ACCESS" HAMPER

- 7/8" high-quality tubular steel with a chrome finish
- Includes one Easy Access bag
- Foot pedal
- 4" non-marking, quiet-ride casters (two locking)
- WARRANTY: Limited

500-HM-1	\$219.99
4" casters	

SINGLE MEDIUM-DUTY HAMPER WITH FOOT PEDAL

- 7/8" high-quality tubular steel with a chrome finish
- Bag lock secures liners in place
- Foot pedal
- 3" non-marking, quiet-ride casters (two locking)
- WARRANTY: Limited

Item #	Price
500-HM-2	\$81.99
3" casters	

#### 3-BUSHEL POLY LAUNDRY HAMPER

- Poly construction with sturdy zinc-plated base
- Leakproof and virtually indestructible
- Foot pedal
- 2" hooded ball casters for smooth transporting
- WARRANTY: Limited

 Item #
 Price

 500-HM-3
 \$157.99

 3-bushel hamper
 \$157.99

#### HAMPER LIDS

Item #	Price
500-HML White	\$26.99
500-HMLID-BLU Blue	\$26.99
500-HMLID-BLK Black	\$26.99
500-HMLID-RED Red	\$26.99
500-HMLID-YEL Yellow	\$26.99

#### HAMPER LABELS

ltem #	Price
500-HMLAB-BLU Blue, "Soiled Linen" lettering	\$36.99
500-HMLAB-WHT White, "Soiled Linen" lettering	\$36.99
500-HMLAB-RED Red, "Infectious Linen" lettering	\$36.99
500-HMLAB-BLK Black, "Biohazardous Waste" lett	\$36.99 ering
500-HMLAB-GRN Green, "Resident Personal Linen"	\$36.99 lettering

#### HAMPER

ltem #	Price
500-HMBAG Nylon	\$38.99
500-HMBAG-A Antimicrobial	\$46.99
500-HMBAG-POLY-YEL Yellow poly disposable infectious li	\$92.99 nen
500-HMBAG-POLY-CLR Clear poly disposable	\$84.99
500-HMBAG-POLY-WHT White poly disposable soiled linen	\$95.99

#### THREE PANEL PRIVACY SCREEN

- 67" overall height
- Each panel measures 27" x 55"
- Three flame-retardant white vinyl panels
- White powder-coated steel frame
- WARRANTY: Limited

ltem #	Price
500-3PPS-A	\$96.95
Privacy screen	
500-3PPS-C	\$114.95
Privacy screen with 2" casters	

# CHECK OUT promedsupply.com

for full descriptions, color options and accessories

# Clinical & Diagnostic

#### WELCH ALLYN<sup>®</sup> <u>SPOT VITAL</u> SIGNS<sup>®</sup>

- Measurements are quick, easy and accurate
- Lightweight and portable
- Large LCD display with backlightIncludes 2 cuffs
- (1 adult size and 1 large adult size)
- Powered by a rechargeable 7.2 volt battery (included)
- WARRANTY: Battery: One year; Spot vital signs unit: Two years

Item #	Price	
WA42NTB-E1 Vital signs unit	\$2,039.99	
WA4700-60 Mobile stand	\$285.99	
WADS44-11 \$41.99 Aneroid sphygmomanometer with adult one-piece cuff		
WADS44-12 Aneroid sphygmomanomete one-piece cuff (large)	\$43.99 er with adult	




WADS44-11

Welch/Allyn<sup>.</sup>

00

1-08

WA42NTB-E1



### **INNOVATIVE PRODUCT**

TEMPERATURE MODULE

# **CUSTOMER FAVORITE**



### FREE SHIPPING ITEMS SHIP WITHIN 48 HOURS

ADVIEW<sup>™</sup> MODULAR DIAGNOSTIC STATION

- Features individual blood pressure, temperature and pulse oximetry modules
- "Connect-and-go" capability between the blood pressure module and the temperature and pulse oximetery modules
- Automated and manual BP options
- Filac<sup>™</sup> FasTemp<sup>™</sup> provides 10-second oral or axillary temperatures
- SpO<sub>2</sub> module includes industryleading Nellcor compatible technology
- WARRANTY: Two years

Item #	Price
500-ST61 Diagnostic station	\$1,699.99
500-ST154 Mobile stand	\$285.99
Item #	Price
500-ST62 Blood pressure module	\$542.77
500-ST63 Temperature and blood pre	\$956.96 ssure modules
500-ST64 Pulse oximetry and blood p	\$1,264.95 ressure modules

# Clinical & Diagnostic

### **INNOVATIVE PRODUCT**

# A TRULY PORTABLE REAL-TIME ULTRASOUND BLADDER SCANNER

#### PORTASCAN<sup>+</sup> LTC BLADDER SCANNER

- Real-time ultrasound scanning and image capture
- Simple, intuitive software with touch-screen input
- 3.5-5.0 MHz abdominal probe for children and adults
- Automatic bladder volume calculation
- Touch-screen length and width calipers for kidney and other organ measurement
- Large database stores up to 250 images, records and reports
- Built-in thermal printer for record keeping and archiving
- Requires a rechargeable battery (included) or external power supply
- WARRANTY: One year

Item #	Price
<b>600-UBSU</b> Bladder scanner system and st	\$7,295.00 tand
600-UBSU-P Replacement convex 3.5-5.0 a	\$2,458.00 bdominal probe
600-UBSU-BA Rechargable battery	\$489.00
600-UBSU-GEL Ultrasonic Gel	\$12.99
600-UBSU-THP Thermal paper	\$79.99

Simple point-and-freeze operation makes the Portascan+ LTC the bladder scanner of choice for the busy clinical environment. The Portascan+ LTC calculates bladder volume automatically from an ultrasound image. A built-in thermal printer provides a permanent hard copy printout of results for accurate recordkeeping and archiving.





CONNECTIVITY

USB INPUT/OUTPUT



#### BRAUN THERMOSCAN<sup>®</sup> PRO 4000 EAR THERMOMETER

- ExacTemp<sup>™</sup> technology ensures positioning and improved accuracy
   Provides fast results
- Automatic probe cover eject button
- Easy-to-read LCD display
- Requires two AA batteries (included)
- WARRANTY: Three-year limited

ltem #	Price	
WA04000-200	\$189.75	
Thermoscan Pro 4000		
WA05074-800	\$17.27	
Disposable probe covers (200/bx 4bx/cs)		

#### WELCH ALLYN<sup>®</sup> SURETEMP<sup>®</sup> PLUS 690 ELECTRONIC THERMOMETER

- Provides 4-second oral, 10-second rectal and 15-second adult axillary temperatures
- Monitor mode capability for measuring 3-5 minute temperatures
- Used probe covers can be ejected
- automaticallyLarge LCD display
- Requires three AA batteries
- (not included)
- WARRANTY: Instrument: Two-year limited; Probe: One-year limited; Probe well: 90-day limited

ltem #	Price
WA690	\$272.99
SureTemp Plus 690	
WA05031-105	\$11.39
Disposable probe covers (	250/bx 20bx/cs)

#### GENIUS<sup>™</sup> 2 THERMOMETER

Item # 150-303000

Genius 2 150-303030

- Determines body temperature directly from the tympanic membrane
- 1-2 second response time
- Programmable display
- Single-use disposable probe coversRequires three AAA batteries

Disposable probe covers (96/bx 22bx/cs)

(included) • WARRANTY: Three years

- NON-CONTACT INFRARED THERMOMETER
- Measure temperature by pointing thermometer at center of forehead
- Temperature measurement is quick, precise, hygienic and safe
- No probe covers necessary
- Requires four AAA batteries (included)
- WARRANTY: Two years

Item #	Price	
600-TH100	\$169.00	
Non-contact infrared thermometer		

	~	1	7	-		
1.57	( A 10			а	54L	13

Item #	Price
500-BA-AA AA alkaline batteries	\$0.42 (EA)
500-BA-AAA AAA alkaline batteries	\$0.42 (EA)

\$285.00

\$6.30

# Clinical & Diagnostic



# HEALTH O METER SCALES

#### FOLDING WHEELCHAIR DIGITAL SCALES

- 1000 lb/454 kg weight capacity
- Features wheels and one or two folding ramps for easy portability and storage
- Wide platform is ideal for wheelchair or bariatric use
- EMR connectivity via USB
- ¾" LCD display
- WARRANTY: Two-year limited

Item #	Price
600-SWPF-1 1 ramp	\$1,698.95
600-SWPF-2 2 ramps	\$1,795.95

#### CHAIR SCALES

 600 lb/270 kg (600-SCHE) or 440 lb/200 kg (600-SCHM) weight capacity

Health o meter

fession

- Two wheels for ease of movement
- Flip arm and foot rest
- WARRANTY: Two-year limited

#### Digital Only

- EMR connectivity via USB
- 1¼″ LCD display

ltem #	Price
600-SCHE	\$499.00
Digital	
600-SCHM	\$429.95
Mechanical	

#### DIGITAL WHEELCHAIR RAMP SCALES

- 1000 lb/454 kg weight capacity
- Dual ramps
- EMR connectivity via USB
- Two wheels for ease of movement
- WARRANTY: Two-year limited

#### PHYSICIAN'S BEAM SCALE

600-SWPP-2

- 390 lb/180 kg weight capacity
- Large, easy-to-read measuring bar
- Built-in height rod
- WARRANTY: 10-year limited

Item #	Price
600-SWPP-1	\$1,399.95
Platform: 31.5" x 31.5" x 2"	
600-SWPP-2	\$1,606.49
Platform: 43" x 42" x 2.5"	

ltem #	Price
600-SPM	\$183.00
390 lb capacity	



#### PHYSICIAN'S BEAM SCALE

- 400 lb/180 kg weight capacity
- Dual-reading die-cast weighbeam
- Heavy-duty platform
- Durable steel construction with an electrostatic powder paint finish
- WARRANTY: One-year limited

#### Item # 500-439 400 lb capacity

\$207.75

#### CHAIR SCALES

- 400 lb/180 kg weight capacity
- Lift-away arms and footrest
- Integral wheels for easy transportRequires six AA batteries (not
- included) or optional AC adapter
- WARRANTY: One-year limited

#### Digital Only

• EMR/EHR ready

Ifem #	Price
<b>500-6475</b> Digital	\$968.85
500-475 Mechanical	\$484.42

#### DIGITAL BARIATRIC PLATFORM SCALE

- 600 lb/270 kg weight capacity
- Low-profile, non-skid platform
- Stainless steel handrails with padded grips
- EMR/EHR ready
- Integral wheels for easy transport
- WARRANTY: One-year limited

**500-6855** \$1,782.61 600 lb capacity

#### PORTABLE WHEELCHAIR SCALES

- 1000 lb/450 kg weight capacity
- For use with wheelchairs, straightlegged chairs or as a stand-on scale
- Sloped ramps for dual-side entry
- EMR/EHR ready
- Requires six C batteries (not included) or optional AC adapter
- WARRANTY: One-year limited

ltem #	Price
<b>500-BRW1000</b> Platform: 40" x 30"	\$2,417.30
<b>500-6550</b> Platform: 32" x 30"; Folding	\$2,083.04

# INTRODUCING A BREAKTHROUGH IN PT/INR TESTING

Coag-Sense PT/INR is the only point of care device that uses **direct detection** of clot formation to test rather than an electrochemical process.

The Coag-Sense PT/INR System provides accurate prothrombin time (PT) results in less than one minute, allowing abnormal results to be addressed immediately – no more losing valuable time waiting for lab results! Its portable tabletop design provides stability during testing and tests can be conducted using fresh capillary whole blood, venous or re-calcified plasma samples.





### FREE METER WITH 5-BOX STARTER KIT

#### COAG-SENSE PT/INR MONITORING SYSTEM

- CLIA waived
- Direct clot detection technology for superb accuracy and precision
- Requires just 10 µ of blood
- Accurate and precise (CV 2.5%) results in one minute
- Measures the actual time required for clotting
- Barcoded test strips eliminate manual calibration
- Requires four AA batteries (included) or an A/C adapter (included)
- WARRANTY: One year with immediate replacement

ltem #	Price
600-PTINR-KIT	\$1,125.00
New Customer Starter Kit includes: Meter Kit (600-PTINR-MTR) and 5 boxes of test strips (600-PTINR50)	
600-PTINR-MTR	\$1,445.00
Meter Kit includes: 1 meter, 1 user manual, 1 quick reference of auto single-usage lancets (quantity: 12), batteries, $1 A/C$ power adapter, 1 minipip 1 rack of pipette tips (quantity: 96)	4 AA alkaline
600-PTINR-HL	\$65.00
5 high and low control strips, 2 control activation solution vials, pack of 14 sample transfer tubes	
600-PTINR-TT PT/INR blood collection tubes (54/pk)	\$10.50
600-PTINR50 PT/INR monitoring test strips (50/bx)	\$289.00
040-970121 Haemolance Plus safety lancets (100/bx 4	\$22.63 .8bx/cs)



- 1 meter (value: \$1,445.00)
- 5 boxes of test strips (quantity: 250; value: \$1,445.00)
- 10 high and low control strips
- 12 auto single-use lancets
- 1 minipipette two-pack
- 96 pipette tips







500-EKG1 & 500-EKGS-2



500-EKG3-R



The CoaguChek XS PT/INR system gives healthcare professionals the flexibility and control they need to confidently manage a wide range of residents on anticoagulation medication. It's designed to provide a high level of quality, accuracy and resident safety.



#### COAGUCHEK XS PT/INR SYSTEM

- CLIA waived
- Accurate results in about one minute
- Requires a small sample size (8 µ)
- Records 300 test results with date and time
- Requires four AAA batteries (included)
- WARRANTY: One-year limited

ltem #	Price
<b>025-15160</b> Test strips with code chip (48/b)	\$245.00
025-5001	\$73.53
CoaquChek lancets (200/bx)	

# FREE METERS WITH PURCHASE

When you commit to bringing at least **five CoaguChek units** into your facility and purchasing **12 boxes of strips for each unit**, you will **receive the meters for free**! For example, if you commit to five meters, you would only need to purchase 60 boxes of strips over the course of up to four years. This gives you the flexibility to purchase strips as you need them rather than paying for all of them up front.

#### BURDICK 8300 ECG/EKG DEVICE

- User-friendly interface
- 5.7" monochrome display
- Manual or automatic 3-, 4-, 6- and 12-lead printouts
- Standard with wired Ethernet and compatible with wireless 802.11
- Requires a 14.4V rechargeable battery pack (included) or AC adapter (included)
- WARRANTY: Parts and labor: Three years

Item #	Price
500-EKG1 ECG/EKG device	\$3,295.00
500-EKGS-2 Stand	\$435.00

# ECG/EKG

#### BURDICK 4250 DIGITAL HOLTER RECORDER

- Streamlined five-button keypad
- Event button allows residents to record episodes that concern them
- 5- or 7-lead 3-channel Holter recordings with an ECG signal preview screen
- Integrated pacemaker detection
- Requires one AAA battery (not included)
- WARRANTY: Unit: One year; Flash card memory: 90 days

Item #	Price
500-EKG3-R	\$1,595.00
Digital holter recorder	

# Respiratory

005-AS098-101 LIMITED TIME PRICE REDUCTION!



CHECK OUT promedsupply.com

for full descriptions on these and other items

#### VISIONAIRE OXYGEN CONCENTRATOR

- Compact and lightweight with near-silent operation
- Weight: 30 lbs
- Oxygen concentration: 1-5 LPM at 90% +5.5/-3%
- 20.8" H x 14.1" W x 11.5" D
- Operational temperature range: 41° F to 105° F
- WARRANTY: *Parts:* Five years

Price
\$473.95
\$507.95

NEWLIFE INTENSITY 10 CONCENTRATOR

- Provides maximum pressure and flow
- Weight: 58 lbs
- Oxygen concentration: 2–9 LPM at 92% ±3%; 10LPM at 90% ±3%

in stati

AIR EP.

- 27.5" H × 16.5" W × 14.5" D
- Operational temperature range: 41° F to 105° F
- WARRANTY: Five years

 Item #
 Price

 005-AS099-10
 \$1,250.00

 NewLife Intensity 10 concentrator

#### FREESTYLE BY AIRSEP

• The only truly wearable portable oxygen concentrator

005-AS099-10

- Weight: 4.4 lbs
- Oxygen concentration: 1-3 pulse settings; equivalent to continuous flow of 90% oxygen ±3%
- 8.6" H x 6.1" W x 3.6" D
- Operational temperature range: 41° F to 104° F
- WARRANTY: Unit: Three years; AirBelt: One year

 Item #
 Price

 005-AS095-1
 \$3,395.00

 Freestyle concentrator

#### 005-AS095-1

500-ASP1

drive

#### HEAVY DUTY PORTABLE SUCTION MACHINE

- Controlled vacuum from 23" HG
- E-Z grip handle
- Base with cord wrap
- UL approved
- WARRANTY: Two-year limited

ltem #	Price
500-ASP1	\$148.95
Suction machine	

#### **CUSTOMER FAVORITE**



- Two-way display
- Not for continuous monitoring
- Rubber boot included
- Requires two AAA batteries (included)
- WARRANTY: Two-year limited

Item #	Price
500-PULSEOXF	\$49.95
Clip style pulse oximeter	
500-BA-AAA	\$0.42 (EA)
AAA alkaline batteries	



- Weight: 4.1 lbs
- Ultra compact, portable design
- Powerful piston pump is efficient, dependable and quiet
- Complete with disposable neb kit
- WARRANTY: Five-year limited

Item #	Price
500-NC01M Pacifica Elite nebulizer	\$19.99





- Customizable audible and visual alarm settings
- Individualized pulse rate and oxygen saturation high/low thresholds
- Durable rubber boot provides protection against accidental drops
- Requires two AA batteries (included)
- WARRANTY: Two-year limited

ltem #	Price
500-PULSEOXH-DRV Includes finger probe	\$249.95
500-BA-AA AA alkaline batteries	\$0.42 (EA)

COMPARTMENT STYLE COMPRESSOR/NEBULIZER



- Maximum pressure: 30 psi
- Operating pressure: 12-14 psi
- Liter flow: 8 L/min
- 45° angle delivery allows resident to recline during treatment
- WARRANTY: Five-year limited

Item # 500-NC01 \$39.95 Compressor/nebulizer

# DEPENDABLE ENTERAL CARE AT YOUR FINGERTIPS



The Compliance Enteral Pump is easy to use and can deliver all enteral feeding formulas. It's simple to program and many of its functions can be activated with a single touch. Reliable and accurate, the pump is designed to give the user access to its many useful features.

#### Additional features:

- Programming directions are located directly on the pump
- Built-in pole clamp for resident mobility
- Convenient one-hand loading
- Large, easy-to-read display windows
- Dose Check records the accumulated amount of formula delivered over several feedings
- Easy-to-differentiate audible and visible alarms
- Memory is preserved in the event of a power failure and when the unit is turned off
- Dose Limit can be set in 1 mL increments up to 2000 mL, then in 5 mL increments up to 9999 mL
- A special cap holder keeps caps close by and helps prevent contamination
- Adjustable display brightness allows for simple day and night operation

#### Ask your ProMed territory manager about our exclusive enteral feeding in-service guide and staff training!

Item #	Description
600-199235C	Enteral Delivery System Feed Pump Compliance OST (1/cs)
600-199001	IV Poles for Compliance Pump Deluxe w/2 Prongs 5 Leg (2/cs)

© Professional Medical, Inc. We reserve the right to correct any errors that may occur within this brochure.

Compliance, Smart Choice, Professional Medical, Inc., ProMed, and We Make Improving Care Easier! are trademarks of Professional Medical, Inc.

Promotions valid only for items listed in catalog at special listed pricing.

Orders will be shipped in the order they are received. Special priced products available while supplies last. Promotional Pricing is valid through September 31, 2013, unless otherwise noted. Orders placed after the promotional expiration date will be sold at our everyday low price. Applicable sales tax will be applied. Our shipping/handling policy procedures apply.

700-Equipment and Furniture Solutions\_Q3\_13 [06/13]



Professional Medical, Inc. A Tradition of Quality, Value, & Trusted Service Since 1968

1917 Garnet Ct. New Lenox, IL 60451 = 800.648.5190 = fax 866.726.7416 = promedsupply.com

We make improving care easier!®